

*Call 1-800-Endoscope!*

**1800Endoscope.com**



## EQUINE ENDOSCOPES – FINDING THE RIGHT ONE.

*When It Comes to Equine Endoscopes, Size Really Does Matter!*

### Tackling the Needs of Today's Equine Vet

Not all horses are created equally and not all equine endoscopes are either. For some procedures, a standard human colonoscope or gastroscope is the perfect fit yet others require a longer reach. In both of these cases, 1800Endoscope.com can help.

### Selection Second to None

Thanks to a climate-controlled on-site warehouse and skilled inventory management, our inventory includes all types of endoscopes. Gastrosopes and colonoscopes are just a few of the types we carry. By maintaining these inventory levels, 1800Endoscope.com offers veterinarians and facilities more selection than most companies. 99% of our inventory is in-stock and immediately available. In addition to the pre-owned equipment, we also carry custom-built 3 meter video and fiber endoscopes. 3 meter scopes are available in several different sizes: **12.8mm O.D.** w/ 3.7mm channel, **11.5mm O.D.** W/ 3.2mm channel, **9.8mm O.D.** W/ 2.8mm channel. These endoscopes can be delivered in two to four weeks.

### New 1800EVFS6135



This New 6mm! diameter fiberscope has 4-way articulation and 135 mm insertion length. 2.0mm Biospy Channel. This scope is used for inspection on medium sized animals, such as dogs, to diagnose blockages. Also for perfect for Equine to examine the guttural pouch. This scope is shipped complete with a detachable ACMI type light

Container, ( This scope uses Syringe to clean lens and insufflate )

### Sales without Service? Not Here

There are two indicators of good service, repeat buyers and referrals. 1800Endoscope.com takes pride in the fact that over 60% of our sales are to repeat buyers and referrals from those same repeat buyers.

### Quality, not Quantity

While 1800Endoscope.com sells a large quantity of endoscopes to the veterinary market, we still believe that quality is the most important factor in maintaining a business. We sell original, pre-owned Olympus®, Pentax®, Fujinon®, and our custom-built brand. All of our endoscopes undergo thorough rigorous testing

### Video Technology 3 Meter Scopes



With one of our 3 Meter Scope system you can perform a multitude of procedures. A few examples are:

**Tieback Surgery**  
**Upper Respiratory Track Surgery**  
**Epiglottic Entrapment**  
**Guttural Pouch Tympany Surgery**  
**Sinuscopy**  
**Laproscopy**

Visit our website at [www.1800Endoscope.com](http://www.1800Endoscope.com). We have a solution for your needs regardless of size, budget, or purpose.

#### Headquarters

6220 Manatee Ave W Suite 301  
Bradenton, FL 34209

[www.1800ENDOSCOPE.com](http://www.1800ENDOSCOPE.com)

Phone: 1-800-363-6726 or 941-792-7138

Fax: 941-761-4613

Email: [sales@1800endoscope.com](mailto:sales@1800endoscope.com)

#### Satellite Office

5757 Blue Lagoon Dr Suite 235  
Miami, FL 33126