

Pre-Owned Olympus®, Pentax®, Fujinon®, and Storz® Specifications

We sell New and Used Endoscopes and are not associated or affiliated in any way with Olympus® Inc., Pentax® Inc., Fujinon® Inc., Machida®, Wilson-Cook® Inc®. All product names and images on our webpages are trademarks or registered ® trademarks of their respective holders, Including 1800Endoscope® LLC, Olympus® Inc., Pentax® Inc., Fujinon® Inc., Wilson-Cook® Inc®. Specifications Subject to change without notice.

For best use, click the Bookmarks Tab

To Search, Use Ctrl-F

Angioscopes

Olympus

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
AF Series											
AF-5	.5 mm	-	60	Forward		150 cm	-	-	-	-	-
AF-8	.8 mm	-	55	Forward		130 cm	-	-	-	-	-
AF-14	1.4 mm	-	75	Forward		130 cm	-	-	-	-	-
????	AF-22A	2.2 mm	-			110 cm					
	AF-28C	2.2 mm	1.0 mm	75	Forward	100 cm		-	-	-	-
PF Series											
PF-8	.8 mm	-	60	Forward		120 cm		-	-	-	-
PF-14	1.4 mm	-	75	Forward		120 cm		-	-	-	-
PF-22	2.2 mm	-	75	Forward		120 cm		-	-	-	-
PF-28	2.7 mm	1.0 mm	75	Forward		120 cm		-	-	-	-

1800Endoscopy.com LLC

Bronchoscopes

Olympus

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
2 Series											
BF-B2											
BF-4B2											
BF-PE2	5.0 mm	2.2 mm	120	Forward	3-50 mm	60 cm	89 cm	180	130	-	-
BF-TE2	6.0 mm	2.8 mm	120	Forward	3-50 mm	60 cm	89 cm	180	130	-	-
10 Series											
BF-10	6.0 mm	2.0 mm	100	Forward	3-50 mm	55 cm	76 cm	160	130	-	-
BF-1T10	6.0 mm	2.6 mm	100	Forward	3-50 mm	55 cm	76 cm	180	100	-	-
BF-2T10	6.0 mm	1.5/2.0 mm	90	Forward	3-50 mm	55 cm	76 cm	160	100	-	-
BF-3C10											
BF-6C10											
BF-P10	5.0 mm	2.0 mm	90	Forward	3-50 mm	55 cm	76 cm	180	100	-	-
20 Series											
BF-1T20	5.9 mm	2.8 mm	100	Forward		55 cm		180	100	-	-
BF-1T20D	6.0 mm	2.8 mm	100	Forward	3-50 mm	55 cm	84 cm	180	100	-	-
BF-20	6.0 mm	2.2 mm	120	Forward		55 cm		180	130	-	-
BF-20D	6.0 mm	2.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
BF-3C20	3.6 mm	1.2 mm	90	Forward		55 cm		160	100	-	-
BF-6C20											
BF-N20	2.2 mm	-	75	Forward	2-50 mm	55 cm	84 cm	160	90	-	-
BF-P20	4.9 mm	2.2 mm	100	Forward		55 cm		180	100	-	-
BF-P20D	5.0 mm	2.2 mm	100	Forward	3-50 mm	55 cm	84 cm	180	100	-	-
BF-XT20	6.4 mm	3.2 mm									
30 Series											
BF-1T30	6.0 mm	2.8 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
BF-30	6.0 mm	2.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
BF-3C30	3.6 mm	1.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-

	BF-P30	5.0 mm	2.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
	BF-XT30	6.3 mm	3.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
	BF-LT30											
40 Series												
	BF-1T40	6.0 mm	2.8 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
	BF-3C40	3.6 mm	1.2 mm	120	Forward	3-50 mm	55 cm	76 cm	180	130	-	-
	BF-40	5.9 mm	2.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
	BF-P40	5.0 mm	2.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
	BF-XP40	2.8 mm	1.2 mm	90	Forward	2-50 mm	60 cm		180	130	-	-
	BF-XT40	6.2 mm	3.2 mm	120	Forward	3-50 mm	55 cm	84 cm	180	130	-	-
60 Series												
	BF-MP60	4.0 mm	2.0 mm	120	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
	BF-P60	4.9 mm	2.2 mm	120	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
	BF-XP60	2.8 mm	1.2 mm	90	Forward	2-50 mm	60 cm		180	130	-	-
	BF-XT60	6.2 mm	3.2 mm	120	Forward		55 cm		180	130	-	-
	BF-1T60	6.0 mm	3.0 mm	120	Forward		60 cm		180	130	-	-
160 Series												
	BF-160	5.2 mm	2.0 mm	120	Forward	3-100 mm	60 cm	87 cm	180	130	-	-
	BF-1T160	6.0 mm	2.8 mm	120	Forward	3-100 mm	60 cm	87 cm	180	130	-	-
	BF-3C160	3.8 mm	1.2 mm	120	Forward	3-100 mm	60 cm	87 cm	180	130	-	-
	BF-MP160F	4.4 mm	2.0 mm	120	Forward	3-50 mm	60 cm	87 cm	180	130	-	-
	BF-P160	4.9 mm	2.0 mm	120	Forward	3-100 mm	60 cm	87 cm	180	130	-	-
	BF-XP160F	2.8 mm	1.2 mm	120	Forward	2-50 mm	60 cm	87 cm	180	130	-	-
	BF-XT160	6.3 mm	3.2 mm	120	Forward	3-100 mm	60 cm	87 cm	180	130	-	-
180 Series												
	BF-Q180-AC	5.3 mm	2.0 mm	120	Forward		60 cm		180	130	-	-
200 Series												
	BF-1T200	6.0 mm	2.6 mm	100	Forward		55 cm		180	130	-	-
	BF-200	5.7 mm	2.0 mm	100	Forward		55 cm		180	130	-	-
	BF-P200	5.2 mm	2.0 mm	100	Forward	3-50 mm	55 cm	82 cm	180	130	-	-

240 Series

BF-1T240	6.0 mm	2.6 mm	120	Forward	3-100 mm	55 cm	82 cm	180	130	-	-
BF-240	5.7 mm	2.0 mm	120	Forward	3-100 mm	55 cm	82 cm	180	130	-	-
BF-6C240	6.0 mm		120	Forward	2-100 mm	55 cm	81 cm	180	130	-	-
BF-P240	5.2 mm	2.0 mm	120	Forward	3-100 mm	55 cm	82 cm	180	130	-	-

260 Series

BF-1T260	6.0 mm	2.8 mm	120	Forward	3-100 mm	60 cm		180	130	-	-
BF-260	4.9 mm	2.0 mm	120	Forward	3-100 mm	60 cm		180	130	-	-
BF-6C260	5.7 mm	2.0 mm	120	Forward		60 cm		180	130	-	-
BF-P260F	4.4 mm	2.0 mm		Forward		60 cm		180	130	-	-
BF-XP260F	2.8 mm	1.2 mm		Forward		60 cm		180	130	-	-

Pentax**BS Series**

FB-15BS	4.8 mm	2.0 mm	100	Forward		60 cm		180	130	-	-
FB-18BS	5.9 mm	2.6 mm	100	Forward		60 cm		180	90	-	-

FB Series

FB-8V	2.7 mm	1.2 mm	100	Forward		60 cm		180	130	-	-
FB-10V	3.4 mm	1.2 mm	95	Forward		60 cm		180	130	-	-
FB-10X	3.5 mm	1.2 mm	95	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
FB-15H											
FB-15P	4.9 mm	2.2 mm	100	Forward		60 cm					
FB-15V	4.9 mm	2.2 mm	120	Forward		60 cm		180	130	-	-
FB-15X	4.9 mm	2.2 mm	100	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
FB-18P	6.0 mm	2.8 mm	100	Forward		60 cm					
FB-18V	5.9 mm	2.8 mm	120	Forward		60 cm		180	130	-	-
FB-18X	6.0 mm	2.2 mm	100	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
FB-18RX	6.0 mm	2.2 mm	120	Forward	3-50 mm	60 cm	90 cm	180	130	-	-
FB-19H											
FB-19TV	6.2 mm	3.2 mm	120	Forward		60 cm		180	130	-	-
FB-19TX	6.2 mm	3.2 mm	100	Forward	3-50 mm	60 cm	90 cm	180	130	-	-

30 Series

EB-1530T												
EB-1530T2												
EB-1530T3	5.3 mm	2.0 mm	120	Forward		60 cm		180	130	-	-	
EB-1830T3	6.0 mm	2.6 mm	120	Forward		60 cm		180	130	-	-	
VB-1530T3												
VB-1830T3												

70 Series

EB-1570	5.1 mm	2.0 mm	120	Forward		60 cm					
EB-1970	6.2 mm	2.8 mm	120	Forward		60 cm					

70K Series

EB-1570K	5.1 mm	2.0 mm	120	Forward		60 cm		210	130	-	-
EB-1970K	6.2 mm	2.8 mm	120	Forward		60 cm		180	130	-	-

Fujinon**BRO Series**

BRO-Y3S	5.7 mm	2.0 mm	100	Forward	3-50 mm	57.5 cm	77 cm	160	100	-	-
BRO-YP2	4.8 mm	2.0 mm	90	Forward	2-50 mm	57.5 cm	77 cm	160	100	-	-
BRO-YL2	6.4 mm	2.6 mm	100	Forward	2-50 mm	57.5 cm	77 cm	150	60	-	-

250 Series

EB-250S	5.3 mm	2.0 mm	120	Forward	3-50 mm	57.5 cm	83 cm	180	130	-	-
---------	--------	--------	-----	---------	---------	---------	-------	-----	-----	---	---

270 Series

EB-270P	3.8 mm	1.2 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-
EB-270S	4.9 mm	2.0 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-
EB-270T	5.9 mm	2.8 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-

310 Series

EB-310S

450 Series

EB-450S	5.3 mm	2.0 mm	120	Forward	3-50 mm	57.5 cm	83 cm	180	130	-	-
EB-450T	6.3 mm	2.8 mm	120	Forward	3-50 mm	57.5 cm	83 cm	180	130	-	-

470 Series

EB-470P	3.8 mm	1.2 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-
EB-470S	4.9 mm	2.0 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-
EB-470T	5.9 mm	2.8 mm	120	Forward	3-100 mm	60 cm	90 cm	180	130	-	-

1800Endoscope.com LLC

Choledochoscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
-------	--------	---------	-----	-----	-----	----	----	----	------	-------	------

Olympus

10 Series

CHF-10											
CHF-P10	5.6 mm	2.0 mm				30 cm					

20 Series

CHF-B20	4.5 mm	1.7 mm	100	Forward	3-50 mm	187 cm	219 cm	160	100	-	-
CHF-CB20	4.6 mm	1.2 mm				70 cm					
CHF-XP20	3.7 mm	1.2 mm	120	Forward	3-50 mm	45 cm	74 cm	160	130	-	-
CHF-XT20											
CHF-P20Q	5.2 mm	2.0 mm	120	Forward	3-50 mm	38 cm	67 cm	160	130	90	90
CHF-P20	4.9 mm	2.2 mm	120	Forward	3-50 mm	38 cm	67 cm	160	130	-	-
CHF-T20	6.0 mm	2.6 mm	120	Forward	3-50 mm	38 cm	67 cm	130	100	-	-

30 Series

CHF-BP30	3.4 mm	1.2 mm	90	Forward	1-50 mm	180 cm	219 cm	130	120	-	-
CHF-CB30L	2.8 mm	1.2 mm	75	Forward	2.5-50 mm	70 cm	101 cm	120	120	-	-
CHF-CB30S	2.8 mm	1.2 mm	75	Forward	2.5-50 mm	45 cm	76 cm	120	120	-	-

Pentax

FCN Series

FCN-15H											
FCN-15X	4.8 mm	2.2 mm	125/83	Forward		35 cm		180	130	-	-

FCP Series

FCP-8P	2.8 mm	.75 mm	90/64	Forward		190 cm		-	-	-	-
--------	--------	--------	-------	---------	--	--------	--	---	---	---	---

FCP-9P	3.1 mm	1.2 mm	90/64	Forward		190 cm		-	-	-	-
--------	--------	--------	-------	---------	--	--------	--	---	---	---	---

Fujinon

CHO Series

CHO-SP	4.9 mm	2.0 mm	90	Forward	2-50 mm	40 cm	70 cm	110	110	-	-
--------	--------	--------	----	---------	---------	-------	-------	-----	-----	---	---

1800Endoscope.com LLC

CF-1T20i	13.7 mm	4.2 mm	120	Forward	5-100 mm	133 cm	165 cm	180	180	160	160
CF-1T20L	13.7 mm	4.2 mm	120	Forward	5-100 mm	168 cm	200 cm	180	180	160	160
CF-20Hi						133 cm					
CF-20HL						168 cm					
CF-20i	13.3 mm	3.2 mm	120	Forward	5-100 mm	133 cm	165 cm	180	180	160	160
CF-20L	13.3 mm	3.2 mm	120	Forward	5-100 mm	168 cm	200 cm	180	180	160	160
CF-20M	13.3 mm	3.2 mm	120	Forward	5-100 mm	103 cm	135 cm	180	180	160	160
CF-P20i	12.2 mm	3.2 mm	120	Forward	5-100 mm	133 cm	165 cm	180	180	160	160
CF-P20L	12.2 mm	3.2 mm	120	Forward	5-100 mm	168 cm	200 cm	180	180	160	160
PCF-20	11.3 mm	2.8 mm	120	Forward	3-100 mm	133 cm	165 cm	180	180	160	160

30 Series

CF-30i	12.9 mm	3.2 mm	140	Forward	5-100 mm	132.5 cm	165 cm	180	180	160	160
CF-30L	12.9 mm	3.2 mm	140	Forward	5-100 mm	167.5 cm	200 cm	180	180	160	160
CF-30M	12.9 mm	3.2 mm	140	Forward	5-100 mm	102.5 cm	135 cm	180	180	160	160

40 Series

CF-40i	12.9 mm	3.2 mm	140	Forward	5-100 mm	133 cm	166.5 cm	180	180	160	160
CF-40L	12.9 mm	3.2 mm	140	Forward	5-100 mm	168 cm	201.5 cm	180	180	160	160

100 Series

CF-100Hi											
CF-100HL											
CF-100i	13.3 mm	3.2 mm	140	Forward	5-100 mm	132 cm	162 cm	180	180	160	160
CF-100L	13.3 mm	3.2 mm	140	Forward	5-100 mm	167 cm	197 cm	180	180	160	160
CF-100Ti	12.9 mm	3.7 mm	140	Forward	5-100 mm	132 cm	162 cm	180	180	160	160
CF-100TL	12.9 mm	3.7 mm	140	Forward	5-100 mm	167 cm	197 cm	180	180	160	160
CF-1T100i	13.7 mm	4.2 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
CF-1T100L	13.7 mm	4.2 mm	140	Forward	5-100 mm	168 cm	198 cm	180	180	160	160
PCF-100	11.3 mm	2.8 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160

130 Series

CF-130i	13.3 mm	3.2 mm	140	Forward		133 cm		180	180	160	160
CF-130L	13.3 mm	3.2 mm	140	Forward		164 cm		180	180	160	160
PCF-130i	11.3 mm	2.8 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
PCF-130L	11.3 mm	2.8 mm	140	Forward	5-100 mm	164 cm	198 cm	180	180	160	160

140 Series

CF-140i	12.9 mm	3.7 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
CF-140L	12.9 mm	3.7 mm	140	Forward	5-100 mm	168 cm	198 cm	180	180	160	160
CF-1T140i	13.7 mm	4.2 mm	140	Forward	5-100 mm	133 cm	198 cm	180	180	160	160
CF-1T140L	13.7 mm	4.2 mm	140	Forward	5-100 mm	168 cm	198 cm	180	180	160	160
CF-2T140											
CF-Q140i	12.9 mm	3.7 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
CF-Q140L	12.9 mm	3.7 mm	140	Forward	5-100 mm	168 cm	198 cm	180	180	160	160
PCF-140i	11.3 mm	3.2 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
PCF-140L	11.3 mm	3.2 mm	140	Forward	5-100 mm	168 cm	198 cm	180	180	160	160

145 Series

CF-Q145i	12.8 mm	3.7 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-Q145L	12.8 mm	3.7 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160

160 Series

CF-2T160i	13.7 mm	3.7/3.2 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-2T160L	13.7 mm	3.7/3.2 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160
CF-Q160Ai	12.8 mm	3.7 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-Q160AL	12.8 mm	3.7 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160
CF-Q160Di	13.2 mm	3.7 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-Q160DL	13.2 mm	3.7 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160
CF-Q160i	12.8 mm	3.7 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-Q160L	12.8 mm	3.7 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160
CF-Q160Zi	12.8 mm	3.7 mm	140/50	Forward	7-100/2-3 mm	133 cm	165 cm	180	180	160	160
CF-Q160ZL	12.8 mm	3.7 mm	140/50	Forward	7-100/2-3 mm	168 cm	200 cm	180	180	160	160
PCF-160Ai	11.5 mm	3.2 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
PCF-160AL	11.5 mm	3.2 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160

165 Series

CF-Q165i	12.8 mm	3.7 mm	140	Forward	3-100 mm	133 cm	165 cm	180	180	160	160
CF-Q165L	12.8 mm	3.7 mm	140	Forward	3-100 mm	168 cm	200 cm	180	180	160	160

180 Series

CF-H180Ai	12.8 mm	3.7 mm	170	Forward	2-100 mm	133 cm	165.5 cm	180	180	160	160
-----------	---------	--------	-----	---------	----------	--------	----------	-----	-----	-----	-----

CF-H180AL	12.8 mm	3.7 mm	170	Forward	2-100 mm	168 cm	200.5 cm	180	180	160	160
CF-Q180Ai	12.8 mm	3.7 mm	170	Forward	3-100 mm	168 cm	200.5 cm	180	180	160	160
CF-Q180AL	12.8 mm	3.7 mm	170	Forward	3-100 mm	133 cm	165.5 cm	180	180	160	160
PCF-Q180Ai	11.5 mm	3.2 mm	140	Forward	3-100 mm	168 cm	200.5 cm	180	180	160	160
PCF-Q180AL	11.5 mm	3.2 mm	140	Forward	3-100 mm	133 cm	165.5 cm	180	180	160	160

200 Series

CF-1T200i											
CF-1T200L											
CF-200HL											
CF-200Hi											
CF-200L											
CF-200i											
CF-200ZL											
CF-200Zi											
CF-2T200L	14.2 mm	3.2/2.8 mm	140	Forward		138 cm		180	180	160	160
PCF-200											

230 Series

CF-230L
CF-230i
PCF-230

240 Series

CF-240Ai	12.0 mm	3.2 mm	140	Forward		133cm		180	180	160	160
CF-240AL	12.0 mm	3.2 mm	140	Forward		168 cm		180	180	160	160
CF-240Di	12.4 mm	3.2 mm	140	Forward		133 cm		180	180	160	160
CF-240DL	12.4 mm	3.2 mm	140	Forward		168 cm		180	180	160	160
CF-240i	12.0 mm	3.2 mm	140	Forward		133 cm		180	180	160	160
CF-240L	12.0 mm	3.2 mm	140	Forward		168 cm		180	180	160	160
CF-Q240Ai	12.9 mm	3.7 mm	140	Forward		133 cm		180	180	160	160
CF-Q240AL	12.9 mm	3.7 mm	140	Forward		168 cm		180	180	160	160
CF-Q240i	12.9 mm	3.7 mm	140	Forward		133 cm		180	180	160	160
CF-Q240L	12.9 mm	3.7 mm	140	Forward		168 cm		180	180	160	160
CF-Q240Zi	12.9 mm	3.2 mm	140/40	Forward		133 cm		180	180	160	160
CF-Q240ZL	12.9 mm	3.2 mm	140/40	Forward		168 cm		180	180	160	160
PCF-240Ai	10.5 mm	3.2 mm	140	Forward		138 cm		180	180	160	160

FC-38SH
 FC-38SX2
 FC-38TLH
 FC-P34
 FC-PP34L

00/01 Series

EC-3400F	11.5 mm	3.5 mm	120	Forward	5-100 mm	150 cm	182 cm	180	180	160	160
EC-3400L	11.5 mm	3.5 mm	120	Forward	5-100 mm	170 cm	202 cm	180	180	160	160
EC-3400M	11.5 mm	3.5 mm	120	Forward	5-100 mm	130 cm	162 cm	180	180	160	160
EC-3800F											
EC-3800L											
EC-3800TL	12.8 mm	3.8/2.8 mm									
EC-3801F	12.8 mm	3.8 mm	120	Forward	5-100 mm	150 cm	182 cm	180	180	160	160
EC-3801F2	12.8 mm	3.8 mm	120	Forward	5-100 mm	150 cm	182 cm	180	180	160	160
EC-3801L	12.8 mm	3.8 mm	120	Forward	5-100 mm	170 cm	202 cm	180	180	160	160
EC-3801L30											
EC-3801M	12.8 mm	3.8 mm	120	Forward	5-100 mm	130 cm	162 cm	180	180	160	160
EC-3801M2	12.8 mm	3.8 mm	120	Forward	5-100 mm	130 cm	162 cm	180	180	160	160
EC-3801S2	12.8 mm	3.8 mm	120	Forward	5-100 mm	115 cm	147 cm	180	180	160	160

30/31 Series

EC-2930F											
EC-3430F	11.7 mm	3.5 mm	120	Forward		150 cm		180	180	160	160
EC-3430L	11.7 mm	3.5 mm	120	Forward		170 cm		180	180	160	160
EC-3430LZ											
EC-3431L											
EC-3830F	12.8 mm	3.8 mm	120	Forward	6-100 mm	150 cm	170 cm	180	180	160	160
EC-3830F2	12.8 mm	3.8 mm	120	Forward	6-100 mm	150 cm	170 cm	180	180	160	160
EC-3830L	12.8 mm	3.8 mm	120	Forward	6-100 mm	170 cm	202.7 cm	180	180	160	160
EC-3830LZ											
EC-3830M	12.8 mm	3.8 mm	120	Forward	6-100 mm	130 cm	162.7 cm	180	180	160	160
EC-3830M2	12.8 mm	3.8 mm	120	Forward	6-100 mm	130 cm	162.7 cm	180	180	160	160
EC-3830S2											
EC-3830TL	12.8 mm	3.8/2.8 mm	120	Forward	6-100 mm	170 cm	203.6 cm	180	180	160	160
EC-3831F	12.8 mm	3.8 mm	140	Forward		150 cm		180	180	160	160
EC-3831L	12.8 mm	3.8 mm	140	Forward		170 cm		180	180	160	160

EC-3831M
EC-3831TL
XEC-3830L

30K/31K Series

EC-2930FK										
EC-3430FK	11.6 mm	3.8 mm	140	Forward	150 cm	180	180	160	160	
EC-3430LK	11.6 mm	3.8 mm	140	Forward	170 cm	180	180	160	160	
EC-3830FK2										
EC-3830LK										
EC-3830MK										
EC-3830MK2										
EC-3830TLK	12.8 mm	3.8/2.8 mm	140	Forward	170 cm	180	180	160	160	
EC-3831FK	12.8 mm	4.2 mm	140	Forward	150 cm	180	180	160	160	
EC-3831LK	12.8 mm	4.2 mm	140	Forward	170 cm	180	180	160	160	
EC-3830LZK										

40 Series

EC-3440F	11.7 mm	3.5 mm	120	Forward	150 cm	180	180	160	160
EC-3440L	11.7 mm	3.5 mm	120	Forward	170 cm	180	180	160	160
EC-3840F	12.8 mm	3.8 mm	120	Forward	150 cm	180	180	160	160
EC-3840F2									
EC-3840L	12.8 mm	3.8 mm	120	Forward	170 cm	180	180	160	160
EC-3840TL	12.8 mm	3.8/2.8 mm	120	Forward	170 cm	180	180	160	160

40K Series

EC-3840LK									
EC-3841LK									
EC-3840TFK	12.8 mm	3.8/2.8 mm	140	Forward	150 cm	180	180	160	160

70K/72K Series

EC-3470LK	11.6 mm	3.8 mm	140	Forward	170 cm	180	180	160	160
EC-3870F2K	12.8 mm	4.2 mm	140	Forward	150 cm	180	180	160	160
EC-3870FK	12.8 mm	4.2 mm	140	Forward	150 cm	180	180	160	160
EC-3870LK	12.8 mm	4.2 mm	140	Forward	170 cm	180	180	160	160
EC-3870LZK	12.8 mm	3.8 mm	140	Forward	170 cm	180	180	160	160
EC-3870M2K	12.8 mm	4.2 mm	140	Forward	130 cm	180	180	160	160

EC-3870MK	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3872LK	12.8 mm	4.2 mm	140	Forward		170 cm		180	180	160	160
EC-3872TLK	12.8 mm	3.8/2.8 mm	140	Forward		170 cm		180	180	160	160

80K Series

EC-3480LK	11.6 mm	3.8 mm	140	Forward		170 cm		180	180	160	160
EC-3880MK	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3880M2K	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3880FK	12.8 mm	4.2 mm	140	Forward		150 cm		180	180	160	160
EC-3880F2K	12.8 mm	4.2 mm	140	Forward		150 cm		180	180	160	160
EC-3880LK	12.8 mm	4.2 mm	140	Forward		170 cm		180	180	160	160

85 Series

EC-3885M	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3885M2	12.8 mm	3.8 mm	120	Forward		130 cm		180	180	160	160
EC-3885F2	12.8 mm	4.2 mm	140	Forward		150 cm		180	180	160	160
EC-3885L	12.8 mm	4.2 mm	140	Forward		170 cm		180	180	160	160

85K Series

EC-3885MK	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3885MK2	12.8 mm	4.2 mm	140	Forward		130 cm		180	180	160	160
EC-3885FK	12.8 mm	4.2 mm	140	Forward		150 cm		180	180	160	160
EC-3885FK2	12.8 mm	4.2 mm	140	Forward		150 cm		180	180	160	160
EC-3885LK	12.8 mm	4.2 mm	140	Forward		170 cm		180	180	160	160
EC-3885TFK	12.8 mm	3.8/2.8 mm	140	Forward		150 cm		180	180	160	160
EC-3885TLK	12.8 mm	3.8/2.8 mm	140	Forward		170 cm		180	180	160	160

Fujinon

Coloscopes

Z Series

FC-1Z	13.0 mm	3.2 mm	120	Forward	3-100 mm	166 cm	198 cm	180	180	160	160
-------	---------	--------	-----	---------	----------	--------	--------	-----	-----	-----	-----

7 Series

EC7-HL2	14.0 mm	3.2 mm	120	Forward	5-100 mm	177 cm	210 cm	180	180	160	160
EC7-HM2	14.0 mm	3.2 mm	120	Forward	5-100 mm	152 cm	185 cm	180	180	160	160

EC7-LR2
EC7-MR2

100 Series

FC-100LR	13.0 mm	3.2 mm	125	Forward	5-100 mm	166 cm	198 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	--------	--------	-----	-----	-----	-----

200 Series

EC-200LR	13.0 mm	3.2 mm	120	Forward	5-100 mm	177 cm	204 cm	180	180	160	160
EC-200LT	13.3 mm	3.7 mm	120	Forward	5-100 mm	177 cm	204 cm	180	180	160	160
EC-200MP(13)	11.1 mm	2.8 mm	120	Forward	5-100 mm	133 cm	160 cm	180	180	160	160
EC-200MR	13.0 mm	3.2 mm	120	Forward	5-100 mm	133 cm	160 cm	180	180	160	160
EC-200MR(15)	13.0 mm	3.2 mm	120	Forward	5-100 mm	152 cm	179 cm	180	180	160	160
EC-200MT(15)	13.3 mm	3.7 mm	120	Forward	5-100 mm	152 cm	179 cm	180	180	160	160
EC-200WM2(15)	13.0 mm	3.2 mm	140	Forward	6-100 mm	152 cm	179 cm	180	180	160	160
EC-200WM2(13)	13.0 mm	3.2 mm	140	Forward	6-100 mm	133 cm	160 cm	180	180	160	160

201 Series

EC-201WI	12.8 mm	3.8 mm	140	Forward	4-100 mm	152 cm	182 cm	180	180	160	160
EC-201WL	12.8 mm	3.8 mm	140	Forward	4-100 mm	169 cm	199 cm	180	180	160	160
EC-201WM	12.8 mm	3.8 mm	140	Forward	4-100 mm	133 cm	163 cm	180	180	160	160

250 Series

EC-250DL5	12.8 mm	3.8/2.8 mm	140	Forward	5-100 mm	169 cm	199 cm	180	180	160	160
EC-250DM5	12.8 mm	3.8/2.8 mm	140	Forward	5-100 mm	133 cm	163 cm	180	180	160	160
EC-250HL2	12.9 mm	3.8 mm	140	Forward	7-100 mm	169 cm	196 cm	180	180	160	160
EC-250HL5	12.8 mm	3.8 mm	140	Forward	5-100 mm	169 cm	199 cm	180	180	160	160
EC-250LP	11.1 mm	3.2 mm	140	Forward	7-100 mm	169 cm	196 cm	180	180	160	160
EC-250LP5	11.1 mm	3.2 mm	140	Forward	5-100 mm	169 cm	199 cm	180	180	160	160
EC-250MP5	11.1 mm	3.2 mm	140	Forward	3-100 mm	133 cm	163 cm	180	180	160	160
EC-250WI5	12.8 mm	3.8 mm	140	Forward	3-100 mm	152 cm	182 cm	180	180	160	160
EC-250WL5	12.8 mm	3.8 mm	140	Forward	3-100 mm	169 cm	199 cm	180	180	160	160
EC-250WM5	12.8 mm	3.8 mm	140	Forward	3-100 mm	133 cm	163 cm	180	180	160	160

300 Series

EC-300CM	14.0 mm	3.2 mm	120	Forward	4-50 mm	133 cm	164 cm	180	180	160	160
EC-300ER	13.0 mm	3.2 mm	120	Forward	5-100 mm	79 cm	110 cm	180	180	160	160
EC-300HL	14.0 mm	3.2 mm	120	Forward	5-100 mm	177 cm	208 cm	180	180	160	160

EC-300HM(15)	14.0 mm	3.2 mm	120	Forward	5-100 mm	152 cm	183 cm	180	180	160	160
EC-300HM(13)	14.0 mm	3.2 mm	120	Forward	5-100 mm	133 cm	164 cm	180	180	160	160
EC-300LR	13.0 mm	3.2 mm	125	Forward	5-100 mm	177 cm	208 cm	180	180	160	160
EC-300MR(15)	13.0 mm	3.2 mm	125	Forward	5-100 mm	152 cm	183 cm	180	180	160	160
EC-300MR(13)	13.0 mm	3.2 mm	125	Forward	5-100 mm	133 cm	164 cm	180	180	160	160
EC-300MS	11.5 mm	2.8 mm	120	Forward	5-100 mm	152 cm	183 cm	180	180	160	160

310 Series

EC-310HL	13.0 mm	3.2 mm	120	Forward	8-100 mm	177 cm	208 cm	180	180	160	160
EC-310HM	13.0 mm	3.2 mm	120	Forward	8-100 mm	133 cm	164 cm	180	180	160	160
EC-310LT	13.0 mm	3.7 mm	120	Forward	3-100 mm	177 cm	208 cm	180	180	160	160

410 Series

EC-410D	14.0 mm	3.2/2.8 mm	140	Forward	6-100 mm	151 cm	182 cm	180	180	160	160
EC-410CL	13.3 mm	3.2 mm	110	Forward	5-20 mm	177 cm	207 cm	180	180	160	160
EC-410CT											
EC-410MP											

450 Series

EC-450WM5	12.8 mm	3.8 mm	140	Forward	3-100 mm	133 cm	163 cm	180	180	160	160
EC-450WI5	12.8 mm	3.8 mm	140	Forward	3-100 mm	152 cm	182 cm	180	180	160	160
EC-450WL5	12.8 mm	3.8 mm	140	Forward	3-100 mm	169 cm	199 cm	180	180	160	160
EC-450HL5	12.8 mm	3.8 mm	140	Forward	3-100 mm	169 cm	199 cm	180	180	160	160
EC-450MP5	11.1 mm	3.2 mm	140	Forward	3-100 mm	133 cm	163 cm	180	180	160	160
EC-450LP5	11.1 mm	3.2 mm	140	Forward	3-100 mm	169 cm	199 cm	180	180	160	160
EC-450DM5	12.8 mm	3.8/2.8 mm	140	Forward	3-100 mm	133 cm	163 cm	180	180	160	160
EC-450DL5	12.8 mm	3.8/2.8 mm	140	Forward	3-100 mm	169 cm	199 cm	180	180	160	160
EC-450ZW5-M	12.8 mm	3.8 mm	140/65	Forward	1.5-3/6-100 mm	133 cm	163 cm	180	180	160	160
EC-450ZW5-L	12.8 mm	3.8 mm	140/65	Forward	1.5-3/6-100 mm	169 cm	199 cm	180	180	160	160

485 Series

EC-485ZW	14.0 mm	3.2 mm	140	Forward	3-6/7-100 mm	133 cm	163.5 cm	180	180	160	160
----------	---------	--------	-----	---------	--------------	--------	----------	-----	-----	-----	-----

490 Series

EC-490WR5	9.3 mm	2.8 mm	140	Forward	6-100 mm	110 cm	140 cm	210	90	100	100
EC-490ZW5-M	12.8 mm	3.8 mm	140/55	Forward	2-3/6-100 mm	133 cm	163 cm	180	180	160	160
EC-490ZW5-L	12.8 mm	3.8 mm	140/55	Forward	2-3/6-100 mm	169 cm	199 cm	180	180	160	160

1800Endoscope.com LLC

Cystoscopes

Olympus

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
CYF Series											
CYF	5.1 mm	2.0 mm	90	Forward	3-50 mm	35.5 cm	66 cm	210	90	-	-
CYF-2	5.3 mm	2.2 mm									
CYF-3	5.2 mm	2.2 mm	120	Forward	3-50 mm	38 cm	67 cm	210	90	-	-
CYF-4	5.2 mm	2.4 mm	120	Forward	3-50 mm	38 cm		210	120	-	-
CYF-4A	5.2 mm	2.4 mm	120	Forward	3-50 mm	38 cm		210	120	-	-
CYF-5	5.5 mm	2.4 mm	120	Forward				210	120	-	-
CYF-5A	5.5 mm	2.4 mm	120	Forward				210	120	-	-
CYF-200											
CYF-240	5.3 mm	2.0 mm	120	Forward				210	120	-	-
CYF-240A	5.3 mm	2.0 mm	120	Forward				210	120	-	-
CYF-V	5.4 mm	2.2 mm	120	Forward				210	120	-	-
CYF-VA	5.4 mm	2.2 mm	120	Forward				210	120	-	-

A = Has suction function

Pentax

FCY Series

FCY-15P2	4.8 mm	2.2 mm	120/80	Forward		40 cm		220	90	-	-
FCY-15RBS	4.9 mm	2.2 mm									

30 Series

ECY-1530	5.3 mm	2.0 mm	125/83	Forward		40 cm		220	120	-	-
----------	--------	--------	--------	---------	--	-------	--	-----	-----	---	---

70 Series

	ECY-1570	5.1 mm	2.0 mm	120/80	Forward		40 cm		220	120	-	-
70K Series	ECY-1570K	5.1 mm	2.0 mm	120/80	Forward		40 cm		220	120	-	-

Fujinon

Cystoscopes

Laprosopes

410 Series

	EL2-TF410	10.5 mm	-	80	Forward	14-120 mm	31 cm	51.5 cm	100	100	90	90
--	-----------	---------	---	----	---------	-----------	-------	---------	-----	-----	----	----

1800Endoscope.com LLC

Duodenoscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
Olympus											
Duodenoscope											
E Series											
JF-E	11.0 mm	2.8 mm	80	Side	5-60 mm	123.5 cm	155.5 cm	120	90	110	90
V Series											
JF-V2	11.0 mm	3.2 mm	100	Side	5-60 mm	123.5 cm	154.5 cm	110	90	110	90
10 Series											
JF-10	10.5 mm	2.2 mm									
JF-1T10	11.0 mm	2.8 mm	80	Side	5-60 mm	123.5 cm	155.5 cm	120	90	110	90
TJF-10	12.5 mm	4.2 mm	80	Side	5-60 mm	124 cm	155.5 cm	120	90	110	90
JF-V10											
TJF-V10											
20 Series											
JF-20	10.5 mm	2.2 mm	80	Side	5-60 mm	124 cm	157 cm	130	120	110	90
JF-M20	13.0 mm	5.5 mm	80	Side	5-60 mm	122 cm	155 cm	120	90	110	90
JF-1T20	11.0 mm	3.2 mm	80	Side	5-60 mm	123.5 cm	155.5 cm	120	90	110	90
TJF-20	12.5 mm	4.2 mm	80	Side	5-60 mm	124 cm	155.5 cm	120	90	110	90
JF-M20											
30 Series											
JF-1T30	11.0 mm	3.2 mm	80	Side	5-60 mm	124 cm	156 cm	120	90	110	90
TJF-30	12.5 mm	4.2 mm	80	Side	5-60 mm	124 cm	156 cm	120	90	110	90
40 Series											
JF-1T40	11.0 mm	3.2 mm	80	Side	5-60 mm	124 cm	156 cm	120	90	110	90
100 Series											

JF-100	11.5 mm	3.2 mm	100	Side	5-60 mm	124.5 cm	154.5 cm	120	90	110	90
TJF-100	12.5 mm	4.2 mm	100	Side	5-60 mm	124.5 cm	154.5 cm	120	90	110	90
130 Series											
JF-130	11.0 mm	3.2 mm	100	Side	5-60 mm	124 cm	154 cm	120	90	110	90
TJF-130	12.5 mm	4.2 mm	100	Side	5-60 mm	124 cm	155 cm	120	90	110	90
140 Series											
JF-140F	11.0 mm	3.2 mm	100	Side	5-60 mm	123.5 cm	153 cm	120	90	110	90
JF-140R	11.0 mm	3.2 mm	100	Side	5-60 mm	124 cm		120	90	110	90
TJF-140F	12.5 mm	4.2 mm	100	Side	5-60 mm	124 cm	154.5 cm	120	90	110	90
TJF-140R											
145 Series											
TJF-145	12.5 mm	4.2 mm	100	Side	5-60 mm	124 cm	155 cm	120	90	110	90
160 Series											
TJF-160F	11.3 mm	4.2 mm	100	Side	5-60 mm	124 cm	155 cm	120	90	110	90
TJF-160R	11.3 mm	4.2 mm	100	Side	5-60 mm	124 cm	155 cm	120	90	110	90
PJF-160	7.5 mm	2.0 mm	100	Side	2-60 mm	124 cm	155 cm	120	90	90	90
TJF-160VF	11.3 mm	4.2 mm	100	Side	5-60 mm	124 cm	155 cm	120	90	110	90
200 Series											
JF-200											
TJF-200											
230 Series											
JF-230											
240 Series											
JF-240											
TJF-240											
260 Series											
JF-260V	11.3 mm	3.7 mm	100					120	90	110	100
TJF-260V	11.3 mm	4.2 mm	100					120	90	110	100

Pentax

FD Series

FD-34V	11.3 mm	4.2 mm	80	Side		125 cm		120	90	110	90
FD-34W	11.3 mm	4.2 mm	80	Side		125 cm		120	90	110	90

00/01 Series

ED-3200	10.8 mm	2.8 mm	100	Side	5-100 mm	125 cm	157 cm	130	130	120	120
ED-3401	11.3 mm	3.8 mm	100	Side	5-100 mm	125 cm	157 cm	130	130	120	120

10/11 Series

ED-3410	11.3 mm	4.2 mm	100	Side	5-100 mm	125 cm	157 cm	130	130	120	120
---------	---------	--------	-----	------	----------	--------	--------	-----	-----	-----	-----

30/31 Series

ED-2330	7.4 mm	2.2 mm	100	Side		103 cm		120	90	110	90
ED-3230	10.8 mm	3.2 mm	110	Side		125 cm		120	90	110	90
ED-3430	11.3 mm	3.8 mm	110	Side		125 cm		120	90	110	90
ED-3430T	11.3 mm	4.2 mm	110	Side		125 cm		120	90	110	90
ED-3630T	12.1mm	4.8 mm	110	Side		125 cm		120	90	110	90

30K/31K Series

ED-3230K	10.8 mm	3.2 mm	110	Side	4-100 mm	125 cm	157 cm	120	90	110	90
ED-3430K	11.3 mm	3.8 mm	120	Side	4-100 mm	125 cm	157 cm	120	90	110	90
ED-3430TK	11.3 mm	4.2 mm	110	Side	4-100 mm	125 cm	157 cm	120	90	110	90

40 Series

ED-3240	10.8 mm	3.2 mm	110	Side		125 cm		120	90	110	90
ED-3440	11.3 mm	3.8 mm	110	Side		125 cm		120	90	110	90
ED-3440T	11.3 mm	4.2 mm	110	Side		125 cm		120	90	110	90

40K Series

ED-3240K	10.8 mm	3.2 mm	110	Side		125 cm					
----------	---------	--------	-----	------	--	--------	--	--	--	--	--

70K/72K Series

ED-3470TK	11.6 mm	4.2 mm	100	Side		125 cm		120	90	105	90
-----------	---------	--------	-----	------	--	--------	--	-----	----	-----	----

	ED-3670TK	12.1 mm	4.8 mm	100	Side		125 cm	120	90	110	90
80K Series											
	ED-3480TK	11.6 mm	4.2 mm	100	Side		125 cm	120	90	110	90
	ED-3680TK	12.1 mm	4.8 mm	100	Side		125 cm	120	90	110	90
85 Series											
	ED-3485T	11.6 mm	4.2 mm	100	Side		125 cm				

Fujinon

7 Series												
	ED7-XL2	10.7 mm	2.8 mm	80	Side	5-70 mm	125 cm	158 cm	130	120	120	90
	ED7-XT2	12.1 mm	4.5 mm	90	Side	5-100 mm	125 cm	158 cm	130	120	120	90
200 Series												
	ED-200XU	12.1 mm	3.7 mm	90	Side	5-100 mm	125 cm	152 cm	130	90	110	90
250 Series												
	ED-250XL5	10.8 mm	3.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90
	ED-250XT5	11.8 mm	4.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90
	ED-250XL8	10.8 mm	3.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90
	ED-250XT8	11.8 mm	4.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90
310 Series												
	ED-310XU											
450 Series												
	ED-450XL5	10.8 mm	3.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90
	ED-450XT5	11.8 mm	4.2 mm	100	Side	4-60 mm	125 cm	155 cm	130	90	110	90

Enteroscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
Olympus											
SIF-SW	5.0 mm	1.0 mm	120	Forward	5-50 mm	280 cm	300 cm	Forward	Forward	Forward	Forward
SIF-10	11.3 mm	2.8 mm	100	Forward	3-100 mm	101 cm		180	100	160	160
SIF-20L	10.5 mm	2.8 mm	120	Forward		218 cm		180	180	160	160
SIF-100	10.9 mm	2.8 mm	140	Forward	5-100 mm	218 cm	248 cm	180	180	160	160
SIF-100L	11.2 mm	2.8 mm	140	Forward	3-100 mm	218 cm	250 cm	180	180	160	160
SIF-Q140	10.5 mm	2.8 mm	140	Forward	3-100 mm	250 cm	250 cm	180	180	160	160
SIF-Q240	9.8 mm	2.8 mm	140	Forward		200 cm		180	180	160	160

Pentax

00/01 Series

VSB-2900	9.8 mm	2.8 mm	100	Forward	5-100 mm	250 cm	282 cm	130	130	130	130
----------	--------	--------	-----	---------	----------	--------	--------	-----	-----	-----	-----

30/31 Series

VSB-3430	11.7 mm	3.5 mm	120	Forward		220 cm		180	180	160	160
----------	---------	--------	-----	---------	--	--------	--	-----	-----	-----	-----

30K/31K Series

VSB-3430K	11.6 mm	3.8 mm	140	Forward		220 cm		180	180	160	160
-----------	---------	--------	-----	---------	--	--------	--	-----	-----	-----	-----

40 Series

VSB-3440	11.7 mm	3.5 mm	120	Forward		220 cm		180	180	160	160
----------	---------	--------	-----	---------	--	--------	--	-----	-----	-----	-----

40K Series

VSB-3440K	11.6 mm	3.8 mm	140	Forward		220 cm		180	180	160	160
-----------	---------	--------	-----	---------	--	--------	--	-----	-----	-----	-----

Fujinon

200 Series

EN-200WM	10.7 mm	2.8 mm	140	Forward	6-100 mm	230 cm	260 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	--------	--------	-----	-----	-----	-----

410 Series

EN-410WM	10.7 mm	2.8 mm	140	Forward	6-100 mm	230 cm	260 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	--------	--------	-----	-----	-----	-----

1800Endoscope.com LLC

Esophogial

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
-------	--------	---------	-----	-----	-----	----	----	----	------	-------	------

Olympus

ENF Series

ENF-10											
ENF-1T10											
ENF-GP	3.6 mm	-	85	Forward	5-50 mm	30 cm	55 cm	130	130	-	-
ENF-L3	4.3 mm	-	85	Forward	3-50 mm	36.5 cm	59 cm	130	130	-	-
ENF-P2	3.6 mm	-									
ENF-P3	3.6 mm	-	85	Forward	5-50 mm	26 cm	49 cm	130	130	-	-
ENF-P4	3.6 mm	-	85	Forward	5-50 mm	30 cm	53 cm	130	130	-	-
ENF-T3	4.8 mm	2.2 mm	85	Forward	5-50 mm	36.5 cm	59 cm	130	130	-	-
ENF-V	3.9 mm	-	90	Forward		36.5 cm		130	130	-	-
ENF-VT	4.8 mm	2.0 mm	90	Forward		36.5 cm		130	130	-	-
ENF-XP	2.2 mm	-	75	Forward	2.5-50 mm	30 cm	53 cm	130	130	-	-
ENF-200											

Pentax

W Series

FE-15W	4.9 mm	2.0 mm	120	Forward		60 cm		210	120	120	120
--------	--------	--------	-----	---------	--	-------	--	-----	-----	-----	-----

40 Series

EE-1540	5.1 mm	2.0 mm	140	Forward		60 cm		210	120	-	-
---------	--------	--------	-----	---------	--	-------	--	-----	-----	---	---

80K Series

EE-1580K	5.1 mm	2.0 mm	140	Forward		60 cm		210	120	-	-
----------	--------	--------	-----	---------	--	-------	--	-----	-----	---	---

Gastrosopes

Olympus

	Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
E Series	GIF-E	9.8 mm	2.8 mm	100	Forward	3-100 mm	103 cm	134.5 cm	210	90	100	100
	GIF-E3	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	134.5 cm	210	90	100	100
	GIF-XPE	7.9 mm	2.0 mm	100	Forward	3-100 mm	103 cm	134.5 cm	210	90	100	100
P Series	GIF-P3	9.0 mm	2.0 mm	100	Forward	3-100 mm	113.5 cm	131.5 cm	210	90	100	100
T Series	GIF-1T	12.6 mm	3.7 mm	100	Forward	3-100 mm	109.5 cm	132 cm	180	90	100	100
V Series	GIF-V GIF-V2	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	133 cm	210	90	100	100
10 Series	GIF-10 GIF-1T10 GIF-1TV10	11.2 mm	3.7 mm	100	Forward		103 cm		210	90	100	100
	GIF-2T10	12.6 mm	3.7/2.8 mm	100	Forward	3-100 mm	101 cm	134.5 cm	180	90	100	100
	GIF-D10 GIF-K10 GIF-P10 GIF-PV10	9.0 mm	2.0 mm	100	Forward	3-100 mm	102.5 cm	134.5 cm	210	90	100	100
	GIF-Q10 GIF-V10 GIF-V10Z	11.0 mm	2.8 mm	120	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-XK10 GIF-XP10	7.9 mm	2.0 mm	100	Forward	3-100 mm	102.5 cm	134.5 cm	210	90	100	100

	GIF-XQ10 GIF-XV10	9.8 mm	2.8 mm	120	Forward	3-100 mm	102.5 cm	134.5 cm	210	90	100	100
20 Series												
	GIF-1T20 GIF-20	11.2 mm	3.7 mm	100	Forward	3-100 mm	102.5 cm	134.5 cm	210	90	100	100
	GIF-2T20	12.6 mm	3.7/2.8 mm	100	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-K20	11.7 mm	2.8 mm	100	30	3-100 mm	103 cm	135 cm	180	100	100	100
	GIF-P20											
	GIF-PQ20	9.0 mm	2.8 mm	100	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-Q20	11.0 mm	2.8 mm	120	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-XK20	10.5 mm	2.0 mm	100	45	3-100 mm	103 cm	135 cm	180	90	90	90
	GIF-XP20	7.9 mm	2.0 mm	100	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-XQ20	9.8 mm	2.8 mm	120	Forward	3-100 mm	102.5 cm	134.5 cm	210	90	100	100
30 Series												
	GIF-1T30	11.2 mm	3.7 mm	100	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-N30	5.3 mm	2.0 mm	120	Forward	3-50 mm	93 cm	125 cm	210	90	100	100
	GIF-P30	9.0 mm	2.2 mm	120	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-Q30	11.0 mm	2.8 mm	120	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-XP30	7.9 mm	2.0 mm	100	Forward		103 cm		210	90	100	100
	GIF-XQ30	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	135 cm	210	90	100	100
	GIF-XT30	13.2 mm	6.0 mm	100	Forward	3-100 mm	101 cm	135 cm	180	100	100	100
40 Series												
	GIF-Q40	11.0 mm	2.8 mm	120	Forward	3-100 mm	103 cm	135.5 cm	210	90	100	100
	GIF-XQ40	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	137 cm	210	90	100	100
100 Series												
	GIF-100	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	133 cm	210	90	100	100
	GIF-1T100	11.7 mm	3.7 mm	120	Forward	3-100 mm	103 cm	133 cm	210	90	100	100
	GIF-2T100	12.6 mm	3.7/2.8 mm	120	Forward	3-100 mm	104 cm	133 cm	210	90	100	100
130 Series												
	GIF-130	9.8 mm	2.8 mm	120	Forward	3-100 mm	103 cm	133 cm	210	90	100	100
	GIF-1T130	11.0 mm	3.7 mm	120	Forward	3-100 mm	104 cm	133 cm	210	90	100	100

230 Series

GIF-N230
 GIF-P230
 GIF-Q230
 GIF-XQ230
 GIF-XP230

240 Series

GIF-Q240
 GIF-Q240Z 9.8 mm 2.8 mm 140/75 Forward 210 90 100 100
 GIF-1T240 11.0 mm 3.7 mm 140 Forward 210 90 100 100
 GIF-2T240 11.8 mm 3.7/2.8 mm 120 Forward 210 90 100 100
 GIF-SP240
 GIF-XK240 9.7 mm 2.8 mm 100 45 180 90 90 90
 GIF-XP240
 GIF-XQ240

260 Series

GIF-H260 9.5 mm 140 Forward 210 90 100 100
 GIF-Q260 9.2 mm 2.8 mm 140 Forward 210 90 100 100
 GIF-XQ260 9.0 mm 2.8 mm 140 Forward 210 90 100 100
 GIF-XP260 6.5 mm 2.0 mm Forward 210 90 100 100

Pentax**FG Series**

FG-15W 4.9 mm 2.0 mm 120 Forward 105 cm
 FG-16V 5.3 mm 2.0 mm 125 Forward 105 cm 180 180 160 160
 FG-24V 7.9 mm 2.0 mm 105 Forward 105 cm 210 120 120 120
 FG-24W 7.9 mm 2.2 mm 105 Forward 105 cm 210 120 120 120
 FG-27X
 FG-29H
 FG-29V 9.8 mm 2.8 mm 100 Forward 105 cm 210 120 120 120
 FG-29W 9.8 mm 2.8 mm 100 Forward 105 cm 210 120 120 120

70K/72K Series

EG-1870K	6.0 mm	2.0 mm	140	Forward		105 cm		210	120	120	120
EG-2470K	8.0 mm	2.4 mm	140	Forward		105 cm		210	120	120	120
EG-2770K	9.0 mm	2.8 mm	140	Forward		105 cm		210	120	120	120
EG-2970K	9.8 mm	3.2 mm	140	Forward		105 cm		210	120	120	120
EG-3470K	11.6 mm	3.8 mm	140	Forward		105 cm		210	120	120	120

80K Series

EG-1580K	5.1 mm	2.0 mm	140	Forward		105 cm		210	120	120	120
EG-2980K	9.8 mm	2.8 mm	140	Forward		105 cm		210	120	120	120
EG-3480K	11.6 mm	3.8 mm	140	Forward		105 cm		210	120	120	120

85 Series

EG-2985	9.8 mm	3.2 mm	140	Forward		105 cm					
EG-3485	11.6 mm	3.8 mm	140	Forward		105 cm					

85K Series

EG-2985K	9.8 mm	2.8 mm	140	Forward		136 cm					
EG-3485K	11.6 mm	3.8 mm	140	Forward		105 cm					
EG-3885TK	12.8 mm	3.8/2.8 mm	140	Forward		105 cm					

Fujinon**Z Series**

FG-1ZP	7.9 mm	2.2 mm	105	Forward	3-100 mm	103 cm	130 cm	210	90	100	100
FG-1Z	9.8 mm	2.8 mm	105	Forward	3-100 mm	103 cm	130 cm	210	90	100	100

7 Series

EG7-FP2	9.8 mm	2.8 mm	105	Forward	5-100 mm	110 cm	143 cm	210	90	100	100
---------	--------	--------	-----	---------	----------	--------	--------	-----	----	-----	-----

UGI Series

UGI-F7	10.0 mm	2.8 mm	105	Forward	5-100 mm	101 cm	135 cm	210	90	100	100
UGI-FP7	9.5 mm	2.8 mm	105	Forward	5-100 mm	101 cm	135 cm	210	90	100	100
UGI-PE7	7.8 mm	2.2 mm	105	Forward	5-100 mm	101 cm	135 cm	210	90	100	100

100 Series

FG-100PE	7.8 mm	2.8 mm	105	Forward	5-100 mm	103 cm	135 cm	210	90	100	100
FG-100CT	12.1 mm	3.7 mm	105	Forward	5-100 mm	116 cm	149 cm	180	90	90	90

200 Series

EG-200CT	12.7 mm	4.2 mm	120	Forward	5-100 mm	110 cm	140cm	210	90	100	100
EG-200FP	10.3 mm	2.8 mm	120	Forward	5-100 mm	110 cm	137 cm	190	90	100	100
EG-200FP(S)	9.8 mm	2.8 mm	120	Forward	5-100 mm	110 cm	137 cm	190	90	100	100
EG-200HR	9.8 mm	2.8 mm	120	Forward	8-100 mm	110 cm	137 cm	190	90	100	100

201 Series

EG-201FP	9.8 mm	2.8 mm	120	Forward	4-100 mm	110 cm	140 cm	210	90	100	100
----------	--------	--------	-----	---------	----------	--------	--------	-----	----	-----	-----

250 Series

EG-250CT5	10.8 mm	3.8 mm	140	Forward	5-100 mm	110 cm	140 cm	210	90	100	100
EG-250D	11.7 mm	3.2/2.8 mm	120	Forward	6-100 mm	109 cm	140 cm	210	90	100	100
EG-250D5	11.5 mm	3.8/2.8 mm	140	Forward	5-100 mm	109 cm	139 cm	210	90	100	100
EG-250PE5	8.1 mm	2.2 mm	120	Forward	4-100 mm	110 cm	140 cm	210	90	100	100
EG-250WR5	9.3 mm	2.8 mm	140	Forward	5-100 mm	110 cm	140 cm	210	90	100	100

270 Series

EG-270N	6.0 mm	2.0 mm	120	Forward	5-100 mm	110 cm	140 cm	210	90	100	100
EG-270N5	5.9 mm	2.0 mm	120	Forward	3-100 mm	110 cm	140 cm	210	90	100	100

300 Series

EG-300CT	11.5 mm	3.7 mm	105	Forward	5-100 mm	110 cm	141 cm	210	90	100	100
EG-300FP	9.8 mm	2.8 mm	120	Forward	5-100 mm	110 cm	141 cm	210	90	100	100
EG-300HR	10.7 mm	2.8 mm	120	Forward	8-100 mm	110 cm	141 cm	210	90	100	100

310 Series

EG-310HR	9.8 mm	2.8 mm	120	Forward	8-100 mm	110 cm	141 cm	210	90	100	100
----------	--------	--------	-----	---------	----------	--------	--------	-----	----	-----	-----

410 Series

EG-410CT	12.7 mm	4.2 mm	120	Forward	5-100 mm	110 cm	140cm	210	90	100	100
----------	---------	--------	-----	---------	----------	--------	-------	-----	----	-----	-----

450 Series

EG-450CT5	10.8 mm	3.8 mm	140	Forward	4-100 mm	110 cm	140 cm	210	90	100	100
EG-450D	11.7 mm	3.2/2.8 mm	120	Forward	6-100 mm	109 cm	140 cm	210	90	100	100
EG-450D5	11.5 mm	3.8/2.8 mm	140	Forward	3-100 mm	109 cm	139 cm	210	90	100	100
EG-450HR											
EG-450PE5	8.1 mm	2.2 mm	120	Forward	4-100 mm	110 cm	140 cm	210	90	100	100
EG-450WR5	9.3 mm	2.8 mm	140	Forward	4-100 mm	110 cm	140 cm	210	90	100	100
EG-450ZW5	9.8 mm	2.8 mm	140/65	Forward	1.5-3/6-100 mm	110 cm	140 cm	210	90	100	100

470 Series

EG-470N	6.0 mm	2.0 mm	120	Forward	5-100 mm	110 cm	140 cm	210	90	100	100
EG-470N5	5.9 mm	2.0 mm	120	Forward	2-100 mm	110 cm	140 cm	210	90	100	100

485 Series

EG-485ZH	10.7 mm	2.8 mm	120	Forward	3-6/7-100 mm	110 cm	140.5 cm	210	90	100	100
----------	---------	--------	-----	---------	--------------	--------	----------	-----	----	-----	-----

1800Endoscope.com LLC

Hysteroscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
Olympus											
HYF Series											
HYF-1T	4.9 mm	2.2 mm	120	Forward	2-50 mm	29 cm	59 cm	120	120	-	-
HYF-P											
HYF-XP	3.1 mm	1.2 mm	100	Forward	2-50 mm	24 cm	54 cm	100	100	-	-
Pentax											
FHY Series											
FHY-10RP	3.5 mm	1.2 mm	110/76	Forward		28 cm		130	130	-	-
FHY-15RP	4.9 mm	2.2 mm	125/83	Forward		28 cm		130	130	-	-
70K Series											
EHY-1570K											

Fujinon

1800Endoscope.com LLC

Intubation Scopes

Olympus

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
LF Series											
LF-1	4 mm	1.2 mm	75	Forward	3-50 mm	60 cm	83 cm	120	120	-	-
LF-2	4.0 mm	1.5 mm	90	Forward	3-50 mm	60 cm	83 cm	120	120	-	-
LF-DP	3.1 mm	1.2 mm	90	Forward	3-50 mm	60 cm	85.5 cm	120	120	-	-
LF-GP	4.1 mm	1.5 mm	90	Forward	4-50 mm	60 cm	85.5 cm	120	120	-	-
LF-P	2.2 mm	-	75	Forward	2-50 mm	60 cm	83 cm	120	120	-	-
LF-T	5.2 mm	2.6 mm	120	Forward	3-50 mm	60 cm	83 cm	180	130	-	-
LF-TP	5.2 mm	2.6 mm	90	Forward	3-50 mm	60 cm	85.5 cm	180	130	-	-
LF-V	4.1 mm	1.2 mm	120	Forward		60 cm		120	120	-	-

Pentax

BS Series

FI-7BS	2.4 mm	-	95	Forward		60 cm		130	130	-	-
FI-9BS	3.0 mm	1.2 mm	90	Forward		60 cm		130	130	-	-
FI-10BS	3.4 mm	1.4 mm	90	Forward		60 cm		130	130	-	-
FI-13BS	4.1 mm	1.8 mm	95	Forward		60 cm		160	130	-	-
FI-16BS	5.1 mm	2.6 mm	95	Forward		60 cm		160	130	-	-

P Series

FI-7P	2.4 mm	-	95	Forward		60 cm		130	130	-	-
FI-10P2	3.4 mm	1.4 mm	90	Forward		60 cm		130	130	-	-
FI-13P	4.2 mm	1.8 mm	95	Forward		60 cm		160	130	-	-

Nasoscopes

Olympus

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
ENF Series											
ENF-1T10											
ENF-GP	3.6 mm	-	85	Forward		30 cm		130	130	-	-
ENF-L3	4.3 mm	-	85	Forward	3-50 mm	36.5 cm	59 cm	130	130	-	-
ENF-P2	3.6 mm	-									
ENF-P3	3.6 mm	-	85	Forward	5-50 mm	26 cm	49 cm	130	130	-	-
ENF-P4	3.6 mm	-	85	Forward	5-50 mm	30 cm	53 cm	130	130	-	-
ENF-T3	4.8 mm	2.2 mm	85	Forward	5-50 mm	36.5 cm	59 cm	130	130	-	-
ENF-V	3.9 mm	-	90	Forward		36.5 cm		130	130	-	-
ENF-VT	4.8 mm	2.0 mm	90	Forward		36.5 cm		130	130	-	-
ENF-V2	3.4 mm	-	90	Forward		30 cm		130	130	-	-
ENF-XP	2.2 mm	-	75	Forward	2.5-50 mm	30 cm	53 cm	130	130	-	-

Pentax

BS Series

FNL-10BS	3.5 mm	-	85	Forward		30 cm		130	130	-	-
FNL-10RBS	3.4 mm	-	75	Forward		30 cm		130	130	-	-

P Series

FNL-7RP3	2.4 mm	-	75	Forward		30 cm		130	130	-	-
FNL-10RAP	3.4 mm	1.2 mm	95	Forward		35 cm		130	130	-	-
FNL-10RP3	3.4 mm	-	75	Forward		30 cm		130	130	-	-
FNL-13RAP	4.1 mm	1.2 mm	100	Forward		35 cm		130	130	-	-
FNL-15RP3	4.8 mm	2.2 mm	75	Forward		30 cm		130	130	-	-

30 Series

VNL-1130	3.8 mm	-	85	Forward		30 cm		130	130	-	-
VNL-1530T	5.3 mm	2.0 mm	85	Forward		30 cm		130	130	-	-

70K Series

VNL-1170K	3.7 mm	-	85	Forward	3-50 mm	30 cm	50 cm	130	130	-	-
-----------	--------	---	----	---------	---------	-------	-------	-----	-----	---	---

Fujinon**NAP Series**

NAP-LS	3.7 mm	-	90	Forward	1-50 mm	41 cm	68 cm	160	90	-	-
NAP-PS	3.7 mm	-	90	Forward	1-50 mm	31 cm	58 cm	160	90	-	-

1800Endoscope.com LLC

Sigmoidoscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
Olympus											
OSF Series											
OSF	13.3 mm	3.2 mm									
OSF-2	12.2 mm	3.2 mm	100	Forward	3-100 mm	60 cm	92 cm	180	180	160	160
OSF-2-35	12.2 mm	3.2 mm	100	Forward	3-100 mm	35 cm	67 cm	180	180	160	160
OSF-3	12.2 mm	3.2 mm	120	Forward	3-100 mm	70 cm	104 cm	180	180	160	160
10 Series											
CF-P10S	12.2 mm	3.2 mm	120	Forward	5-100 mm	63 cm	95 cm	180	180	160	160
CF-V10S											
20 Series											
CF-P20S	12.2 mm	3.2 mm	120	Forward	5-100 mm	63 cm	95 cm	180	180	160	160
60 Series											
OSF-V60	11.3 mm	3.2 mm	145	Forward	3-100 mm	73 cm		180	180	160	160
100 Series											
CF-100S	13.3 mm	3.2 mm	140	Forward	5-100 mm	62 cm	92 cm	180	180	160	160
130 Series											
CF-130S	13.3 mm	3.2 mm	140	Forward	5-100 mm	64 cm	93.5 mm	180	180	160	160
140 Series											
CF-140S	13.3 mm	3.7 mm	140	Forward	5-100 mm	73 cm	103 cm	180	180	160	160
160 Series											
CF-Q160S	13.2 mm	3.7 mm	140	Forward	3-100 mm	73 cm	104 cm	180	180	160	160
200 Series											
CF-200S	13.3 mm	3.2 mm	140	Forward		62 cm		180	180	160	160

240 Series

PCF-240S	13.3 mm	3.2 mm	140	Forward		62 cm		180	180	160	160
----------	---------	--------	-----	---------	--	-------	--	-----	-----	-----	-----

260 Series**Pentax****FS Series**

FS-34A	11.5 mm	3.5 mm	95	Forward	3-100 mm	70 cm	87 cm	180	180	100	100
FS-34B	12.8 mm	3.8 mm									
FS-34P	11.5 mm	3.5 mm									
FS-34P2	11.5 mm	3.5 mm									
FS-34V	11.5 mm	3.5 mm	120	Forward		70 cm		180	180	160	160
FS-34W	11.5 mm	3.5 mm	120	Forward		70 cm		180	180	160	160

00/01 Series

ES-3801	12.8 mm	3.8 mm	120	Forward	5-100 mm	70 cm	102 cm	180	180	160	160
---------	---------	--------	-----	---------	----------	-------	--------	-----	-----	-----	-----

30/31 Series

ES-3830	12.8 mm	3.8 mm	120	Forward	6-100 mm	70 cm	102.7 cm	180	180	160	160
ES-3831	12.8 mm	3.8 mm	140	Forward		70 cm		180	180	160	160

30K/31K Series

ES-3831K	12.8 mm	4.2 mm	140	Forward		70 cm		180	180	160	160
----------	---------	--------	-----	---------	--	-------	--	-----	-----	-----	-----

40 Series

ES-3840	12.8 mm	3.8 mm	120	Forward		70 cm		180	180	160	160
---------	---------	--------	-----	---------	--	-------	--	-----	-----	-----	-----

70K/72K Series

ES-3870K	12.8 mm	4.2 mm	140	Forward		70 cm		180	180	160	160
----------	---------	--------	-----	---------	--	-------	--	-----	-----	-----	-----

80K Series

ES-3880K	12.8 mm	4.2 mm	140	Forward	70 cm	180	180	160	160
----------	---------	--------	-----	---------	-------	-----	-----	-----	-----

85K Series

ES-3885K	12.8 mm	4.2 mm	140	Forward	70 cm				
----------	---------	--------	-----	---------	-------	--	--	--	--

Fujinon

100 Series

ES-100ER	13.0 mm	3.2 mm	125	Forward	5-100 mm	69.5 cm	101.5 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	---------	----------	-----	-----	-----	-----

200 Series

ES-200ER	13.0 mm	3.2 mm	140	Forward	6-100 mm	79 cm	109 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	-------	--------	-----	-----	-----	-----

250 Series

ES-250WE5	12.8 mm	3.8 mm	140	Forward	5-100 mm	79 cm	109 cm	180	180	160	160
-----------	---------	--------	-----	---------	----------	-------	--------	-----	-----	-----	-----

300 Series

ES-300ER											
----------	--	--	--	--	--	--	--	--	--	--	--

410 Series

ES-410WE	13.0 mm	3.2 mm	140	Forward	6-100 mm	79 cm	109 cm	180	180	160	160
----------	---------	--------	-----	---------	----------	-------	--------	-----	-----	-----	-----

450 Series

ES-450WE5	12.8 mm	3.8 mm	140	Forward	4-100 mm	79 cm	109 cm	180	180	160	160
-----------	---------	--------	-----	---------	----------	-------	--------	-----	-----	-----	-----

1800Endoscope.com LLC

Uteroscopes

Model	IT Dia	Channel	FOV	DOV	DOF	WL	TL	Up	Down	Right	Left
-------	--------	---------	-----	-----	-----	----	----	----	------	-------	------

Olympus

URF Series

URF											
URF-P	3.6 mm	1.2 mm	75	Forward	3-50 mm	70 cm	101 cm	160	100	-	-
URF-P2	3.3 mm	1.2 mm	90	Forward	1-50 mm	70 cm	101 cm	180	100	-	-
URF-P3	2.8 mm	1.2 mm	90	Forward	2-50 mm	70 cm	101 cm	180	180	-	-

Pentax

FUR Series

FUR-9P	3.0 mm	1.2 mm	90/64	Forward		70 cm		180	130	-	-
FUR-9RBS											
FUR-13P	4.1 mm	1.8 mm	110/76	Forward		70 cm		180	130	-	-

1800Endoscope.com LLC

Rigid - Manufacturer

Manufacturer	Model	Dia	Channel	DOV	DOF	WL	TL	Video
ACMI	FL-8168M	4 mm	-	30				
ACMI	FO-8168D	4 mm	-	0				
ACMI	FO-8168L	4 mm	-	70				
ACMI	LAP 10-0 FDA	10 mm						
ACMI	LAP 10-0FD	10 mm	-					
ACMI	LAP 10-0W	10 mm	-					
ACMI	LAP 10-O FD	10 mm						
ACMI	LAP 10-OA	10 mm						
ACMI	LAP 5-0 M	5 mm	-	0				
ACMI	LAP 5-30 M	5 mm						
ACMI	M2-12	4 mm	-	12				
ACMI	M2-30	4 mm	-	30				
ACMI	M3-30	4 mm	-	30				
ACMI	M3-70	4 mm	-	70				
ACMI	M3-70	4 mm	-	70				
ACMI	M3G-0	4 mm	-	0				
ACMI	M3G-12	4 mm	-	12				
ACMI	M3G-30	4 mm	-	30				
Dyonics	3622	4 mm	-	30				
Dyonics	3623	4 mm	-	70				
Dyonics	3626	4 mm	-	30				Y
Dyonics	3628	4 mm	-	70				Y
Dyonics	3768	2.7 mm	-	30				Y
Dyonics	3797	2.7 mm	-	70				Y
Dyonics	3863	1.9 mm	-	30				Y
Dyonics	3896	4 mm	-	30				
Dyonics	3952	4 mm	-	70				Y
Dyonics	3953	4 mm	-	70				Y
Dyonics	4125	4 mm	-	30				Y
Dyonics	4131	2.7 mm	-	30				Y
Dyonics	4339	4 mm	-	30				
Dyonics	4486	4 mm	-	30		15.7 cm		Y
Dyonics	8000	10 mm	-	0				
Dyonics	8008	10 mm		0				
Dyonics	8009	5 mm	-					
Dyonics	7205827	4 mm	-	30				
Dyonics	7205922	4 mm	-	30				Y
Dyonics	7208327	5 mm	-	30				
Linvatec	7403	4 mm	-	30				
Linvatec	7404	4 mm	-	70				
Linvatec	7431	2.7 mm	-	30				
Linvatec	7447	2.4 mm	-	30				
Linvatec	7500	4 mm	-	30				
Linvatec	7501	4 mm	-	70				
Linvatec	7502	4 mm	-	90				
Linvatec	7503	4 mm	-	30				
Linvatec	7560	4 mm	-	30				Y
Linvatec	7561	4 mm	-	70				

Linvatec	7583	4 mm	-	30	
Linvatec	7403R	4 mm	-	30	
Linvatec	7503R	4 mm	-	30	
Linvatec	7560R	4 mm	-	30	
Linvatec	7583R	4 mm	-	30	
Linvatec	C3200	4 mm	-	30	
Linvatec	T4301	4 mm	-		
Lynvatec	7574	5 mm	-	0	
Olympus	A1932A	4 mm	-		
Olympus	A1933A	4 mm	-		
Olympus	A2010	4 mm	-	0	
Olympus	A2011	4 mm	-		
Olympus	A2011A	4 mm	-		
Olympus	A2012A	4 mm	-	30	
Olympus	A2013	4 mm	-	70	
Olympus	A2013A	4 mm	-	70	
Olympus	A2014A	4 mm	-	110	
Olympus	A2015A	4 mm	-	30	
Olympus	A2031A	4 mm	-	12	
Olympus	A2033A	5 mm	-		
Olympus	A3105	4 mm	-		
Olympus	A3405	5 mm	-		
Olympus	A4605	5 mm	-		
Olympus	A4605A	4 mm	-	30	
Olympus	A4673A	3 mm	-	12	
Olympus	A5212	5 mm	-	5	
Olympus	A5240A	10 mm	-	0	26.5 cm
Olympus	A5254A	10 mm	-	0	
Olympus	A5258A	5 mm	-		
Olympus	A5287A	10 mm	-	0	30 cm
Olympus	A7500A	2.5 mm	-	0	
Olympus	A7541	4 mm	-	30	
Olympus	A7590	4 mm	-	0	14 cm
Olympus	A7591	4 mm	-	30	
Olympus	A7592	4 mm	-	70	
Olympus	A7594A	4 mm	-	0	
Olympus	A75975A	4 mm	-	30	
Olympus	A7597A	4 mm	-	70	
Richard Wolf	8471.42	4 mm	-	25	
Richard Wolf	8471.445	4 mm	-	70	
Richard Wolf	8471.45	4 mm	-	70	
Richard Wolf	8475.42	4 mm	-	25	
Richard Wolf	8475.45	4 mm	-	70	
Richard Wolf	8670.31	2.7 mm	-	170	
Richard Wolf	8851.42	5 mm	-	25	
Richard Wolf	8853.42	5 mm	-	25	17 cm
Richard Wolf	8860.433	1.9 mm	-	25	
Richard Wolf	8880.402	4 mm	-	25	
Richard Wolf	8880.431	4 mm	-		
Richard Wolf	8880.433	4 mm	-	25	
Richard Wolf	8880.533	4 mm	-	25	
Richard Wolf	8880.543	4 mm	-	25	

Richard Wolf	8933.441	7.2 mm	-		
Richard Wolf	8934.40	10 mm	-	0	
Richard Wolf	8934.441	10 mm	-	0	
Richard Wolf	8935.431	4 mm	-	0	
Richard Wolf	8950.31	2.7 mm	-	5	45 cm
Richard Wolf	8950.31	2.7 mm	-		
Richard Wolf	8950.35	2.7 mm	-	70	45 cm
Storz	10320B	5.5 mm	-	30	
Storz	10328BA	4 mm	-	30	
Storz	12015A	5.5 mm	-	0	
Storz	26003AA	10 mm	-	0	
Storz	26003BA	10 mm	-		
Storz	26003FA	10 mm	-		
Storz	26006A	5 mm	-		
Storz	26006AA	5 mm	-	0	
Storz	26033AP	10 mm	-	0	
Storz	26033AR	10 mm	-		
Storz	2603APA	10 mm	-	0	
Storz	27005AA	4 mm	-		
Storz	27005FIA	4 mm	-	12	
Storz	27015A	4 mm	-	0	
Storz	27015C	4 mm	-		
Storz	27018B	2.7 mm	-		
Storz	27020A	2.7 mm	-	0	
Storz	27023AA	2.7 mm	-	6	45 cm
Storz	27023CB	2.7 mm	-	70	45 cm
Storz	27029B	2.4 mm	-	30	
Storz	28720D	4 mm	-	90	
Storz	28721AWA	4 mm	-	0	
Storz	7200A	4 mm	-		
Storz	7200B	4 mm	-	30	
Storz	7200BW	4 mm	-	30	
Storz	7200C	4 mm	-	70	
Storz	7200CW	4 mm	-	70	
Storz	7210AA	4 mm	-	0	
Storz	7210BW	4 mm	-	30	
Storz	7210BWA	4 mm	-	30	
Storz	7210CA-DR	4 mm	-	70	
Storz	7210CKA	4 mm	-	70	
Storz	7230FA	4 mm	-	45	
Storz	8702D	5.5 mm	-	90	
Storz	8710CA	5 mm	-	70	
Storz	8712BA	5 mm	-	30	
Storz	8712CA	5 mm	-	70	
Storz	8714BA	5 mm	-	30	
Stryker	0502-385-030	5 mm	-	30	
Stryker	0502-457-045	10 mm	-	45	
Stryker	0502-585-030	5 mm	-	30	
Stryker	0502-777-030	4 mm	-	30	
Stryker	357-010	10 mm	-	0	
Stryker	380-10	10 mm	-	0	
Stryker	380-11	10 mm	-	0	

Stryker	385-30	5 mm	-	30	
Styker	303-231	4 mm	-	30	Y
Styker	307-031	4 mm	-	30	Y
Styker	327-031	4 mm	-	30	Y
Styker	327-071	4 mm	-	70	Y
Styker	334-30	4 mm	-	30	
Styker	334-70	4 mm	-	70	
Styker	342-030	2.7 mm	-	30	
Styker	342-10	2.4 mm	-	10	
Styker	342-11	2.4 mm	-	0	
Styker	343-031-001	2.7 mm	-	30	
Styker	344-31	4 mm	-	30	
Styker	344-71	4 mm	-	70	
Styker	377-032	4 mm	-	30	
Styker	377-31	4 mm	-	30	
Styker	377-32	4 mm	-	30	
Styker	377-72	4 mm	-	70	
Styker	427-031	4 mm	-	30	
Styker	502-142-030	1.9 mm	-	30	
Styker	502-243-030	2.7 mm	-	30	
Styker	502-343-030	2.7 mm	-	30	
Styker	502-427-303	4 mm	-	30	
Styker	502-477-011	4 mm	-	0	
Styker	502-477-031	4 mm	-	30	
Styker	502-527-030	4 mm	-	30	
Styker	502-543-010	2.7 mm	-	0	
Styker	7-343-231	2.7 mm	-	30	
Styker	7-366-31	4 mm	-	30	
Styker	7-366-71	4 mm	-	70	
Styker	7-377-32	4 mm	-	30	
Styker	7-427-31	4 mm	-	30	Y

18000Endoskope.com

Olympus® Light Source and Processor Compatibility

Type/Series

Model

Angioscope

AF Series

AF-5
AF-8
AF-14
AF-22A
AF-28C

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	OTV-A	CV-100	CV-140	CV-160	CV-200	CV-240
			X		Optional					
			X		Optional					
			X		Optional					
			X		Optional					
			X		Optional					

PF Series

PF-8
PF-14
PF-22
PF-28

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	OTV-S6	CV-100	CV-140	CV-160	CV-200	CV-240
		X		X						
		X		X						
		X		X						
		X		X						

Bronchoscopes

2 Series

BF-B2
BF-4B2
BF-PE2
BF-TE2

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						

10 Series

BF-10
BF-1T10
BF-2T10
BF-3C10
BF-6C10
BF-P10

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						

20 Series

BF-1T20
BF-1T20D
BF-20
BF-20D

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						
X	X	X	X	X						

BF-3C20	X	X	X	X	X						
BF-6C20	X	X	X	X	X						
BF-N20	X	X	X	X	X						
BF-P20	X	X	X	X	X						
BF-P20D	X	X	X	X	X						
BF-XT20	X	X	X	X	X						

30 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
BF-1T30		X	X	X	X	X						
BF-30		X	X	X	X	X						
BF-3C30		X	X	X	X	X						
BF-P30		X	X	X	X	X						
BF-XT30		X	X	X	X	X						
BF-LT30		X	X	X	X	X						

40 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
BF-1T40		X	X	X	X	X						
BF-3C40		X	X	X	X	X						
BF-40		X	X	X	X	X						
BF-P40		X	X	X	X	X						
BF-XP40		X	X	X	X	X						
BF-XT40		X	X	X	X	X						

60 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
BF-MP60		X	X	X	X	X						
BF-P60		X	X	X	X	X						
BF-XP60		X	X	X	X	X						
BF-XT60		X	X	X	X	X						
BF-1T60		X	X	X	X	X						

160 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
BF-160				X				X	X	X		
BF-1T160				X				X	X	X		
BF-3C160				X				X	X	X		
BF-MP160F				X				X	X	X		
BF-P160				X				X	X	X		
BF-XP160F				X				X	X	X		
BF-XT160				X				X	X	X		

20 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-1T20i	X	X	X									
CF-1T20L	X	X	X									
CF-20Hi	X	X	X									
CF-20HL	X	X	X									
CF-20i	X	X	X									
CF-20L	X	X	X									
CF-20M	X	X	X									
CF-P20i	X	X	X									
CF-P20L	X	X	X									
PCF-20	X	X	X									

30 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-30i	X	X	X									
CF-30L	X	X	X									
CF-30M	X	X	X									

40 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-40i	X	X	X									
CF-40L	X	X	X									

100 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-100Hi			X				X	X	X			
CF-100HL			X				X	X	X			
CF-100i			X				X	X	X			
CF-100L			X				X	X	X			
CF-100Ti			X				X	X	X			
CF-100TL			X				X	X	X			
CF-1T100i			X				X	X	X			
CF-1T100L			X				X	X	X			
PCF-100			X				X	X	X			

130 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-130i			X				X	X	X			
CF-130L			X				X	X	X			
PCF-130i			X				X	X	X			
PCF-130L			X				X	X	X			

140 Series

CF-140i
 CF-140L
 CF-1T140i
 CF-1T140L
 CF-2T140
 CF-Q140i
 CF-Q140L
 PCF-140i
 PCF-140L

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
		X				X	X	X		
		X				X	X	X		
		X				X	X	X		
		X				X	X	X		
		X					X	X		
		X					X	X		
		X				X	X	X		
		X				X	X	X		

145 Series

CF-Q145i
 CF-Q145L

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
		X					X	X		
		X					X	X		

160 Series

CF-2T160i
 CF-2T160L
 CF-Q160Ai
 CF-Q160AL
 CF-Q160Di
 CF-Q160DL
 CF-Q160i
 CF-Q160L
 CF-Q160Zi
 CF-Q160ZL
 PCF-160Ai
 PCF-160AL

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
		X				X	X	X		
		X				X	X	X		
		X					X	X		
		X					X	X		
		X					X	X		
		X					X	X		
		X					X	X		
		X					X	X		
		X					X	X		
		X				X	X	X		
		X				X	X	X		

plus CLV-160

165 Series

CF-Q165i
 CF-Q165L

Light Sources					Processors					
CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
		X					X	X		
		X					X	X		

plus CLV-160

plus CV-180

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
140 Series	JF-140F			X			X	X	X			
	JF-140R			X			X	X	X			
	TJF-140F			X			X	X	X			
	TJF-140R			X			X	X	X			
145 Series	TJF-145	Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
				X				X				
160 Series	TJF-160F TJF-160R PJF-160 TJF-160VF	Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
				X				X	X	X		
				X				X	X	X		
		X				X	X	X				
		X				X	X	X				
200 Series	JF-200 TJF-200	Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
				X				X			X	
				X				X			X	
230 Series	JF-230	Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
				X					X		X	
240 Series	JF-240 TJF-240	Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
				X						X	X	
				X						X	X	
260 Series	JF-260V TJF-260V	Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
				X				X	X	X	X	
				X				X	X	X	X	

30 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-1T30	X	X	X								
	GIF-N30	X	X	X								
	GIF-P30	X	X	X								
	GIF-Q30	X	X	X								
	GIF-XP30	X	X	X								
	GIF-XQ30	X	X	X								
	GIF-XT30	X	X	X								

40 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-Q40	X	X	X								
	GIF-XQ40	X	X	X								

100 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-100			X				X	X	X		
	GIF-1T100			X				X	X	X		
	GIF-2T100			X				X	X	X		

130 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-130			X				X	X	X		
	GIF-1T130			X				X	X	X		

140 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-140			X				X	X	X		
	GIF-1T140			X				X	X	X		
	GIF-1TQ140Y			X					X	X		
	GIF-P140			X				X	X	X		
	GIF-Q140			X					X	X		
	GIF-XQ140			X					X	X		

145 Series		Light Sources					Processors					
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
	GIF-Q145			X					X	X		

160 Series		Light Sources			Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200
	GIF-160			X			X	X	X		
	GIF-1T160			X			X	X	X		
	GIF-1TQ160			X				X	X		
	GIF-2T160			X			X	X	X		
	GIF-Q160			X				X	X		
	GIF-Q160Z			X				X	X		
	GIF-XP160			X			X	X	X		
	GIF-XTQ160			X				X	X		

plus CLV-160

165 Series		Light Sources			Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200
	GIF-Q165			X				X	X		

180 Series		Light Sources			Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200
	GIF-H180			X					X		
	GIF-N180			X					X		
	GIF-Q180			X					X		

plus CLV-160

200 Series		Light Sources			Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240
	GIF-200			X				X		X	X
	GIF-200Z			X				X		X	X
	GIF-2T200			X				X		X	X
	GIF-2TK200			X				X		X	X
	GIF-Q200			X				X		X	X
	GIF-XK200			X				X		X	X

230 Series		Light Sources			Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240
	GIF-N230			X					X		X
	GIF-P230			X					X		X
	GIF-Q230			X					X		X
	GIF-XQ230			X					X		X
	GIF-XP230			X					X		X

240 Series		Light Sources				Processors					
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240
	GIF-Q240			X						X	X
	GIF-Q240Z			X						X	X
	GIF-1T240			X				X	X	X	X
	GIF-2T240			X				X	X	X	X
	GIF-SP240			X					X	X	X
	GIF-XK240			X				X	X	X	X
	GIF-XP240			X					X	X	X
	GIF-XQ240			X						X	X

260 Series		Light Sources				Processors					
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240
	GIF-H260										X
	GIF-Q260									X	X
	GIF-XQ260									X	X
	GIF-XP260								X	X	X

Hysteroscopes	HYF Series		Light Sources				Processors						
			CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
		HYF-1T	X	X	X	X	X						
		HYF-P	X	X	X	X	X						
		HYF-XP	X	X	X	X	X						

Intubation Scopes	LF Series		Light Sources				Processors						
			CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	OTV-S7	CV-100	CV-140	CV-160	CV-200	CV-240
		LF-1	X	X	X	X	X						
		LF-2	X	X	X	X	X						
		LF-DP	X	X	X	X	X						
		LF-GP	X	X	X	X	X						
		LF-P	X	X	X	X	X						
		LF-T	X	X	X	X	X						
		LF-TP	X	X	X	X	X						
		LF-V			X			X					

Nasoscopes	ENF Series		Light Sources				Processors						
			CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	OTV-S7	CV-100	CV-140	CV-160	CV-200	CV-240
		ENF-1T10	X	X	X	X	X						
		ENF-GP	X	X	X	X	X						
		ENF-L3	X	X	X	X	X						

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	OTV-S7	CV-100	CV-140	CV-160	CV-200	CV-240
ENF-P2		X	X	X	X	X						
ENF-P3		X	X	X	X	X						
ENF-P4		X	X	X	X	X						
ENF-T3		X	X	X	X	X						
ENF-V				X			X					
ENF-VT				X			X					
ENF-V2				X			X					
ENF-XP		X	X	X	X	X						

Sigmoidoscopes

OSF Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
OSF		X	X	X								
OSF-2		X	X	X								
OSF-2-35		X	X	X								
OSF-3		X	X	X								

10 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-P10S		X	X	X								
CF-V10S		X	X	X								

20 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-P20S		X	X	X								

60 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-60	CV-100	CV-140	CV-160	CV-200	CV-240
OSF-V60				X			X					

100 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-100S				X				X	X	X		

130 Series

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
CF-130S				X				X	X	X		

		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
140 Series	CF-140S			X				X	X	X		
		Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
160 Series	CF-Q160S			X					X	X		
		Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
200 Series	CF-200S			X					X		X	X
		Light Sources				Processors						
		CLK-3/4	CLV	CLV-U20/40	CLV-A	CLV-S20	CV-100	CV-140	CV-200	CV-230	CV-240	CV-260
240 Series	PCF-240S			X					X	X	X	X
Uteroscopes												
		Light Sources				Processors						
		CLK-3/4	CLV	CLV-10/U20/40	CLV-A	CLV-S20	CV-1	CV-100	CV-140	CV-160	CV-200	CV-240
URF Series	URF	X	X	X	X	X						
	URF-P	X	X	X	X	X						
	URF-P2	X	X	X	X	X						
	URF-P3	X	X	X	X	X						

1800Endoscopy.com LLC

Pentax® Light Source and Processor Compatibility

Bronchoscopes

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
BS Series	FB-15BS							X		
	FB-18BS							X		

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
FB Series	FB-8V								X	
	FB-10V								X	
	FB-10X								X	
	FB-15H								X	
	FB-15P								X	
	FB-15V								X	
	FB-15X								X	
	FB-18P								X	
	FB-18V								X	
	FB-18X								X	
	FB-18RX								X	
	FB-19H								X	
	FB-19TV								X	
FB-19TX								X		

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30 Series	EB-1530T			X		X	X			
	EB-1530T2			X		X	X			
	EB-1530T3			X		X	X			
	EB-1830T3			X		X	X			
	VB-1530T3			X		X	X			
	VB-1830T3			X		X	X			

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
FC-38LW							X	
FC-38LX							X	
FC-38MW							X	
FC-38MW2							X	
FC-38MX							X	
FC-38SH							X	
FC-38SX2							X	
FC-38TLH							X	
FC-P34							X	
FC-PP34L							X	

00/01 Series

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
EC-3400F		X	X	X	X			
EC-3400L		X	X	X	X			
EC-3400M		X	X	X	X			
EC-3800F		X	X	X	X			
EC-3800L		X	X	X	X			
EC-3800TL		X	X	X	X			
EC-3801F		X	X	X	X			
EC-3801F2		X	X	X	X			
EC-3801L		X	X	X	X			
EC-3801L30		X	X	X	X			
EC-3801M		X	X	X	X			
EC-3801M2		X	X	X	X			
EC-3801S2		X	X	X	X			

30/31 Series

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
EC-2930F		X		X	X			
EC-3430F		X		X	X			
EC-3430L		X		X	X			
EC-3430LZ		X		X	X			
EC-3431L		X		X	X			
EC-3830F		X		X	X			

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
EC-3830F2		X		X	X			
EC-3830L		X		X	X			
EC-3830LZ		X		X	X			
EC-3830M		X		X	X			
EC-3830M2		X		X	X			
EC-3830S2		X		X	X			
EC-3830TL		X		X	X			
EC-3831F		X		X	X			
EC-3831L		X		X	X			
EC-3831M		X		X	X			
EC-3831TL		X		X	X			
XEC-3830L		X		X	X			

30K/31K Series

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
EC-2930FK	X	X						
EC-3430FK	X	X						
EC-3430LK	X	X						
EC-3830FK2	X	X						
EC-3830LK	X	X						
EC-3830MK	X	X						
EC-3830MK2	X	X						
EC-3830TLK	X	X						
EC-3831FK	X	X						
EC-3831LK	X	X						
EC-3830LZK	X	X						

40 Series

Processor/Light Source Combo						Light Source Only		
EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
EC-3440F		X		X	X			
EC-3440L		X		X	X			
EC-3840F		X		X	X			
EC-3840F2		X		X	X			
EC-3840L		X		X	X			
EC-3840TL		X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
00/01 Series	ED-3200			X	X	X	X			
	ED-3401			X	X	X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
10/11 Series	ED-3410					X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30/31 Series	ED-2330			X		X	X			
	ED-3230			X		X	X			
	ED-3430			X		X	X			
	ED-3430T			X		X	X			
	ED-3630T			X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30K/31K Series	ED-3230K	X	X							
	ED-3430K	X	X							
	ED-3430TK	X	X							

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40 Series	ED-3240			X		X	X			
	ED-3440			X		X	X			
	ED-3440T			X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40K Series	ED-3240K	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
70K/72K Series	ED-3470TK	X	X							
	ED-3670TK	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
80K Series	ED-3480TK	X	X							
	ED-3680TK	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
85 Series	ED-3485T									

Enteroscopes

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
00/01 Series	VSB-2900			X	X	X	X			

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30/31 Series	VSB-3430			X		X	X			

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30K/31K Series	VSB-3430K	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40 Series	VSB-3440			X		X	X			

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40K Series	VSB-3440K	X	X							

Esophogioscopes

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
W Series	FE-15W									

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40 Series	EE-1540									

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
80K Series	EE-1580K									

Gastrosopes

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
FG Series	FG-15W								X	
	FG-16V								X	

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
	FG-24V								X	
	FG-24W								X	
	FG-27X								X	
	FG-29H								X	
	FG-29V								X	
	FG-29W								X	
	FG-34W								X	

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
00/01 Series	EG-2700			X	X	X	X			
	EG-2901			X	X	X	X			
	EG-3400			X	X	X	X			
	EG-3800T			X	X	X	X			
	EG-3800TL			X	X	X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30/31 Series	EG-2430			X		X	X			
	EG-2530			X		X	X			
	EG-2730			X		X	X			
	EG-2731			X		X	X			
	EG-2930			X		X	X			
	EG-2931			X		X	X			
	EG-3431			X		X	X			
	EG-3831T			X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30K/31K Series	EG-2730K	X	X							
	EG-2930K	X	X							
	EG-2931K	X	X							
	EG-3430K	X	X							
	EG-3830TK	X	X							

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40 Series	EG-1540					X	X			
	EG-1840					X	X			
	EG-2540					X	X			
	EG-2940					X	X			
	EG-3440					X	X			
	EG-3840T					X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40K Series	EG-3840TK	X	X							



		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
70K Series	EHY-1570K	X	X							

Intubation

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
BS Series	FI-7BS							X		
	FI-9BS							X		
	FI-10BS							X		
	FI-13BS							X		
	FI-16BS							X		

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
P Series	FI-7P							X	X	
	FI-10P2							X	X	
	FI-13P							X	X	

Nasoscopes

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
BS Series	FNL-10BS							X		
	FNL-10RBS							X		

		Processor/Light Source Combo					Light Source Only			
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
P Series	FNL-7RP3							X	X	
	FNL-10RAP							X	X	
	FNL-10RP3							X	X	
	FNL-13RAP							X	X	
	FNL-15RP3							X	X	

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30 Series	VNL-1130			X		X	X			
	VNL-1530T			X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
70K Series	VNL-1170K	X	X							

Sigmoidoscopes

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
FS Series	FS-34A								X	
	FS-34B								X	
	FS-34P								X	
	FS-34P2								X	
	FS-34V								X	
	FS-34W								X	

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
00/01 Series	ES-3801			X	X	X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30/31 Series	ES-3830			X		X	X			
	ES-3831			X		X	X			

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
30K/31K Series	ES-3831K	X	X							

		Processor/Light Source Combo				Light Source Only				
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
40 Series	ES-3840					X	X			

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
70K/72K Series	ES-3870K	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
80K Series	ES-3880K	X	X							

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
85K Series	ES-3885K	X								

Uteroscopes

		Processor/Light Source Combo						Light Source Only		
		EPK-700	EPK-1000	EPM-1000	EPM-3000	EPM-3300	EPM-3500	BS-LH2	LH-150PC	66111
FUR Series	FUR-9P								X	
	FUR-9RBS								X	
	FUR-13P								X	

1800Endoscopy.com

Fujinon® Light Source and Processor Compatibility

Bronchoscopes

BRO Series

BRO-Y3S
BRO-YP2
BRO-YL2

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
X	X	X	X												
X	X	X	X												
X	X	X	X												

250 Series

EB-250S

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
					X										X

270 Series

EB-270P
EB-270S
EB-270T

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
															X
															X
															X

310 Series

EB-310S

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
								X							

450 Series

EB-450S
EB-450T

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
													X		
													X		

470 Series

EB-470P
EB-470S
EB-470T

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
													X		
													X		
													X		

250 Series

EC-250DL5
 EC-250DM5
 EC-250HL2
 EC-250HL5
 EC-250LP
 EC-250LP5
 EC-250MP5
 EC-250WI5
 EC-250WL5
 EC-250WM5

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X
				X	X										X

300 Series

EC-300CM
 EC-300ER
 EC-300HL
 EC-300HM(15)
 EC-300HM(13)
 EC-300LR
 EC-300MR(15)
 EC-300MR(13)
 EC-300MS

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			
						X	X	X	X	X	X	X			

310 Series

EC-310HL
 EC-310HM
 EC-310LT

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
								X	X	X	X				
								X	X	X	X				
								X	X	X	X				

410 Series

EC-410D
 EC-410CL
 EC-410CT
 EC-410MP

Light Sources				Processors - ALL EPX											
FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
													X	X	
													X	X	
													X	X	
													X	X	

450 Series		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
	EC-450WM5														X	X	
	EC-450WI5														X	X	
	EC-450WL5														X	X	
	EC-450HL5														X	X	
	EC-450MP5														X	X	
	EC-450LP5														X	X	
	EC-450DM5														X	X	
	EC-450DL5														X	X	
	EC-450ZW5-M														X	X	
	EC-450ZW5-L														X	X	

485 Series		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
	EC-485ZW														X		

490 Series		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
	EC-490WR5														X		
	EC-490ZW5-M														X		
	EC-490ZW5-L														X		

Cysto/Laprosopes

410 Series		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
	EL2-TF410														X		

Duodenoscopes

7 Series		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
	ED7-XL2							X	X	X	X	X	X				
	ED7-XT2							X	X	X	X	X	X				

270 Series	EG-270N EG-270N5	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
						X	X										X
						X	X										X

300 Series	EG-300CT EG-300FP EG-300HR	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
								X	X	X	X	X	X				
								X	X	X	X	X	X				
								X	X	X	X	X	X				

310 Series	EG-310HR	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
										X			X				

410 Series	EG-410CT	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
															X		

450 Series	EG-450CT5 EG-450D EG-450D5 EG-450HR EG-450PE5 EG-450WR5 EG-450ZW5	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
															X	X	
															X	X	
															X	X	
															X	X	
															X	X	
															X	X	

470 Series	EG-470N EG-470N5	Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
															X		
															X		

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
485 Series	EG-485ZH														X		

Nasoscopes

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
NAP Series	NAP-LS	X	X	X	X												
	NAP-PS	X	X	X	X												

Sigmoidoscopes

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
100 Series	ES-100ER							X	X	X	X	X	X	X			

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
200 Series	ES-200ER					X	X										X

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
250 Series	ES-250WE5					X	X										X

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
300 Series	ES-300ER							X	X	X	X	X	X	X			

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
410 Series	ES-410WE														X		

		Light Sources				Processors - ALL EPX											
		FIL-150	FIL-150EEB	FLX-300	FH-NS	200	201	301A	301B	301C	302A	302B	302C	304	VP-401	VP-450	2200
450 Series	ES-450WE5														X	X	

Olympus Common Terms/Acronyms

Prefixes		Suffixes		Inside P/N	
AF	Angioscope	2 Numbers	Fiber	A	Suction Function
BF	Bronchoscope			D	
CF	Colonoscope	ex. BF-10/CF-10L/GIF-P20		I	Intermediate (Colono 133cm)
CF(S)	Sigmoidoscope			K	
CHF	Choledochoscope	3 Numbers	Video	L	Long (Colono 168cm)
CLE	Light Source			M	Short (Colono 103cm)
CLH	Light Source	ex. CF-100L/GIF-140		P	Pediatric
CLK	Light Source			Q	Super Image
CLV	Light Source			T	Therapeutic
CV	Processor			V	Video
CYF	Cystoscope			X	Slim
ENF	ENT Scope				
GIF	Gastroscope				
HYF	Hysteroscope				
JF	Duodenscope				
LF	Intubation Scope				
LTF	Rigid - Deflectable Tip				
OSF	Sigmoidoscope				
PCF	Pediatric Colonoscope				
PF	Angioscope				
SIF	Small Bowel Scope				
SU	Suction Pump				
TJF	Therapeutic Duoden				
URF	Ureteroscope				
XLS	Light Source				

1800Endoscope.com LLC

Pentax Common Terms/Acronyms

Prefixes	Inside	Suffixes
Type F = Fiber V = Video E = Electronic Video Application B = Bronchoscope C = Colonoscope CN = Choledoch/Nephroscope CP = Cholangio (Baby) scope CY = Cystoscope D = Duodenoscope E = Esophagoscope G = Gastroscope I = Intubation Scope NL = Nasalaryngoscope S = Sigmoidoscope SB = Small Bowel Enteroscope UR = Ureteroscope	Outer Diameter Size = ex. 38/3 = 12.7 mm Processor 30/31/40 = EPM-1000/3300/3500 30K/31K/70K = EPK-700	Scope Characteristics AP = Air Pulse Technology F = Full (1500mm/150cm) L = Long (1700mm/170cm) Q = Oblique Optics R = Rigid Ocular T = Therapeutic/Two Channel U = Ultrasound - Linear UR = Ultrasound - Radial Z = Zoom Series BS = Portable (Bedside) K = Color CCD (Requires EPK-700) P = Simplified Fiberscopes V = Fiber Endoscope

2=Stiffer

Examples using above info:

EC-3830L	EG-2931K	FI-7P
E = Electronic Video C = Colonoscope 38 = 12.7 mm Dia 30 = EPM-1000/3300/3500 L = Long (1700mm/170cm)	E = Electronic Video G = Gastroscope 29 = 9.7 mm Dia 31K = EPK-700 K = Color CCD (See Above)	F = Fiber I = Intubation Scope 7 = 2.4 mm P = Simplified Series

2/4 numbers

Fujinon Common

M	Short
I	Intermed
L	Long
W	Wide
H	Super Image
D	Dual Channel
T	Ther

1800Endoscope.com LLC