

EG-201FP Specifications

Viewing direction	Forward	
Observation range	4 ~ 100mm	
Field of view	120 °	
Distal end diameter	9.8mm	
Flexible portion diameter	9.8mm	
Bending capability :UP	210 °	
DOWN	90 °	
LEFT	100 °	
RIGHT		
Forceps channel diameter	2.8mm	
Working length	1,100mm	
Total length	1,400mm	



EC-201WM/WI/WL Specifications

Viewing direction	Forward	
Observation range	4 ~ 100mm	
Field of view	140 °	
Distal end diameter	12.8mm	
Flexible portion diameter	12.8mm	
Bending capability :UP	180 °	
DOWN		
LEFT	160 °	
RIGHT		
Forceps channel diameter	3.8mm	
Working length	1,330/1,520/1,690mm	
Total length	1,630/1,820/1,990mm	

EPX-201H Specifications

TV system ——— NTSC/PAL
 Image output ——— Video(BNC) 1 channel
 S-Video 3 channels
 Functions ——— • Electronic magnification X1.25, X1.5(2steps)
 • Image Enhancement • Color bar display
 • Color control • Illuminaiton control
 • Metering mode: Average and Peak measurement
 • Patients data recording for 23 patients
 Power supply ——— AC120V/230V
 Light source ——— 150 watt Halogen lamp
 Dimention ——— 420(W) X 450(D) X 160(H)mm, 16.5 X 17.7 X 6.3 inches
 Weight ——— 25kg (55lbs)



FUJI PHOTO OPTICAL CO., LTD.

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN
 TEL:81-48-668-2153 FAX:81-48-668-1570
 E-mail:sales@msv.fujinon.co.jp URL http://www.fujinon.co.jp/

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
 TEL:1-973-633-5600 FAX:1-973-633-8818
 E-mail:med.info@fujinon.com URL http://www.fujinon.com/

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY
 TEL:49-2154-924-0 FAX:49-2154-924-290
 E-mail:Fujinon@fujinon.de URL http://www.fujinon.de/

FUJI KOKI SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD #10-04, HENDERSON INDUSTRIAL PARK,
 SINGAPORE 159552
 TEL:65-6276-4988 FAX:65-6276-6911 E-mail:fujinon@fujinon.com.sg

F.P.FOCUSING PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W.2128, AUSTRALIA
 TEL:61-2-9748-2744 FAX:61-2-9748-2428 E-mail:fpfocus@dot.net.au

Shanghai FUJINON Medical Equipment Co. Ltd.

FLAT B, 6/FL, SHI MEI BLD, 79 RI JING ROAD,
 WAI GAO QIAO FREE TRADE ZONE, SHANGHAI 200131, CHINA
 TEL:86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)
 FAX:86-21-5424-9264 E-mail:sh-fujinon@etang.com