## 1800Endoscope.com®

## 1800PVS9150

Portable Veterinary Scope



The 1800PVS series is used for visualizing areas where the human eye cannot see. This product uses a mini CCD (Charged Coupled Device). The PVS series has 4 way articulation. There are many great features that come standard with the new PVS. Our new optical system has high resolving power and sensitivity which makes for a good quality image and is adjustable through the user friendly operating system. The light-weight wireless design makes it very easy to handle and clean with its new water proof design.

## **Specifications**

CCD color Image Sensor Resolution 325K x 250K Depth of view:

10~100mm

**Angulations:** 

100 up

90 down

90 left

90 right

Illumination: White LED

**Outer Diameter of** 

insertion tube: 9mm

Working Length: 160cm

Image Output: AV,USB

Power: 3.7V 2200mAh

(lithium battery)

**Audio: Integrated** 

microphone

Biopsy channel: 2.0mm



1800Endoscope.com 1.800.363.6726 941.792.7138 941.761.4613 Fax Sales@1800endoscope.com